IMPORTANT INFORMATION – PLEASE READ BEFORE USE  
INFORMATIONS IMPORTANTES – À LIRE AVANT UTILISATION  
INFORMAZIONI IMPORTANTI – LEGGERE PRIMA DELL'USO  
WICHTIGE INFORMATIONEN – VOR GEBRAUCH BITTE LESEN  
INFORMACIÓN IMPORTANTE – LÉALA ANTES DE UTILIZAR ESTE PRODUCTO  
重要な情報 - 使用前にお読みください  
ВАЖНА ИНФОРМАЦИЯ – МОЛЯ, ПРОЧЕТЕТЕ ПРЕДИ УПОТРЕБА  
VAŽNE INFORMACIJE – PROČITATI PRIJE UPORABE  
DŮLEŽITÉ INFORMACE – PŘEČTĚTE SI PŘED POUŽITÍM  
VIGTIGE OPLYSNINGER – SKAL LÆSES INDEN BRUG  
BELANGRIJKE INFORMATIE – VOOR GEBRUIK DOORLEZEN  
OLULINE TEAVE – PALUN LUGEGE ENNE KASUTAMIST  
TÄRKEITÄ TIETOJA – LUETTAVA ENNEN KÄYTTÖÄ  
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ – ΠΑΡΑΚΑΛΟΥΜΕ ΔΙΑΒΑΣΤΕ ΠΡΙΝ ΑΠΟ ΤΗ ΧΡΗΣΗ  
FONTOS TUDNIVALÓK – HASZNÁLAT ELŐTT OLVASSA EL  
SVARĪGA INFORMĀCIJA — LŪDZU, IZLASIET PIRMS LIETOŠANAS  
SVARBI INFORMACIJA – PERSKAITYKITE PRIEŠ NAUDODAMI  
VIKTIG INFORMASJON – LES FØR BRUK  
WAŽNE INFORMACJE – NALEŻY PRZECZYTAĆ PRZED UŻYCIEM  
INFORMAÇÃOES IMPORTANTES — LER ANTES DE UTILIZAR  
INFORMAȚII IMPORTANTE — VĂ RUGĂM SĂ LE CITIȚI ÎNAINTEA UTILIZĂRII  
DŮLEŽITÉ INFORMÁCIE – PREČÍTAJTE SI ICH PŘED POUŽITÍM  
ПОМЕМБНЕ ИНФОРМАЦИЈЕ – ПРЕБЕРИТЕ ПРЕД УПОРАБО  
VIKTIG INFORMATION – LÄS FÖRE ANVÄNDNING  
ÖNEMLI BİLGİ - LÜTFEN KULLANMADAN ÖNCE OKUYUN

Contents (English - EN).....	1
Contenu (Français - FR).....	6
Contenuto (Italiano - IT) .....	12
Inhalt (Deutsch - DE) .....	18
Índice (Español - ES).....	24
内容 (日本語 - JA).....	30
Съдържание (Български - BG) .....	35
Sadržaj (Hrvatski - HR) .....	41
Obsah (Česky - CS) .....	46
Indhold (Dansk - DA).....	51
Inhoud (Nederlands - NL).....	57
Sisukord (Eesti - ET) .....	63
Sisältö (Suomi - FI) .....	69
Περιεχόμενα (Ελληνικά - EL).....	75
Tartalom (Magyar - HU) .....	81
Saturs (Latviski - LV) .....	87
Turinys (Lietuvių - LT).....	92
Innhold (Norsk - NO) .....	98
Spis treści (Polski - PL) .....	103
Índice (Português - PT).....	109
Cuprins (Română - RO).....	115
Obsah (Slovenčina - SK).....	121
Vsebina (Slovenščina - SL).....	127
Innehåll (Svenska - SV) .....	133
İçindekiler (Türkçe - TR).....	138

## INSTRUCTIONS FOR USE

### **Intended Use**

Tecomet delivery systems (cases and trays) are intended for the protection, organization, and the delivery of surgical instruments and/or other medical devices to the surgical field. Tecomet delivery systems (cases and trays) are not designed to maintain sterility by themselves. They are designed to facilitate the sterilization process when used in conjunction with a wrapping material (FDA cleared sterilization wrap) or a specified filtered sterilization container system. Wrapping materials and sterilization containers are designed to allow air removal, steam penetration/evacuation (drying), and maintain the sterility of the internal components.

### **Recommendations for Care, Cleaning, and Sterilization of Tecomet Delivery Systems (cases and trays)**

Both physical and chemical (detergent) processes may be necessary to clean soiled items. Chemical (detergent) cleaners alone cannot remove all soil and debris; therefore, a careful manual cleaning of each item with a soft sponge or cloth is essential for maximum decontamination. For difficult access areas, a clean soft bristled brush is recommended. Once the items have been cleaned, they should be thoroughly rinsed with clean water to remove any detergent or chemical residue before sterilization. Tecomet recommends the use of a mild enzymatic detergent with a near neutral pH. Do not use solvents, abrasive cleaners, metal brushes, or abrasive pads. Cases and trays may be placed in mechanical cleaning equipment.

**Always** inspect for cleanliness or damage before use. Make sure all latches and handles are secure and in working order. Do not overload cases and trays. Balance contents uniformly within the container and arrange to allow steam to come in contact with all objects in the container. Tecomet recommends that the delivery systems (cases and trays) be processed according to the sterilization wrap or filtered container manufacturer's instructions prior to sterilization to maintain sterility of internal components/items and for proper aseptic presentation to the surgical field.

For steam autoclaves, Tecomet has verified product performance in the following cycles. Our testing was done with common surgical appliances such as reamers, drills, hammers, rasps, drivers, chisels, awls, handpieces, rongeurs, blades and bits, including short lumened devices.

Complex instruments may require disassembly and extended sterilization times. Always follow instrument manufacturer instructions if their sterilization or drying recommendations exceed these guidelines.

### **1. Prevacuum Sterilizer**

Wrapped cases, trays, and instruments or cases, trays, and instruments should be exposed to 132°C to 135° (270°F to 275°F) for at least four (4) minutes. Dry for twenty (20) to forty (40) minutes.

### **2. Gravity Displacement Sterilizer**

Wrapped cases, trays, and instruments should be exposed to 132°C to 135°C (270°F to 275°F) for at least thirty (30) minutes, or 121°C to 123°C (250°F to 254°F) for at least fifty-five (55) minutes. Dry for twenty (20) to fifty (50) minutes.

Variables that may affect drying times include: loading density of the case/tray, instrument configuration, total contents of the sterilizer, steam quality, equipment maintenance and others.

### **3. Ethylene Oxide Sterilization**

Tecomet delivery systems (cases and trays) have been validated for Ethylene Oxide Sterilization using 10/90 cycle (10%EO/90%HCFC). A minimum two hour exposure cycle at 52°C to 57°C (125°F to 135°F) is required. Our validation testing was done with common surgical appliances, scopes, and scope accessories, including stainless steel lumened devices up to 8½ inches. Complex instruments may require disassembly and extended sterilization times. Always follow instrument manufacturer instructions if their sterilization recommendations exceed these guidelines. Always follow instrument manufacturer instructions regarding use in Ethylene Oxide systems. Cases and trays have been shown to be within acceptable residual ranges after degassing for twelve (12) hours.

#### **4. Filtered Sterilization Container Systems**

Tecomet delivery systems (cases and trays) have been validated for placement inside Ultra Container, Genesis, and Aesculap filtered sterilization container systems and exposed to 132°C to 135°C (270°F to 275°F) for at least four (4) minutes for prevacuum sterilizers only. Tecomet delivery systems have not been validated for use in filtered sterilization containers in gravity displacement sterilizers and should not be used in such a system.

Follow the recommendations of the manufacturer of the filtered sterilization containers for the proper placement and use of cases and trays inside the container.

#### **Materials and Restricted Substance**

- For indication that the device contains a restricted substance or material of animal origin see product label.

#### **Contraindications**

Tecomet delivery systems (cases and trays) have **NOT** been validated for use with flexible endoscopes, or devices with lumens or working channels longer than 4 inches (3mm ID). Always refer to instrument manufacturer instructions.

Tecomet delivery systems (cases and trays) have **NOT** been validated for ETO sterilization of devices with lumens or working channels longer than 8½ inches (1.5mm ID). Always refer to instrument manufacturer instructions.

Stacking of delivery systems (cases and trays) and overloading of the units will adversely affect sterilization and drying effectiveness. **DO NOT STACK** cases and trays in autoclave chamber.

#### **Notice:**

1. **DO NOT** load cases into sterilizer on sides or upside down with lid side on the shelf or cart. Load cases on cart or shelf, so that the lid is always facing upward. This will allow for proper drying. Tecomet cases and trays are designed to drain in this position.
2. After the autoclave door is opened, all cases must be allowed to cool thoroughly. Place cases on a rack or shelf with linen cover until cooling is complete. The potential for condensation may increase if the case is not allowed to cool properly.

3. If condensation is observed, check to ensure steps 1 and 2 were followed. In addition, verify that the steam, which is being used for sterilization processing, has a quality of more than 97%. Also confirm that the sterilizers have been inspected for routine maintenance in accordance with manufacturer recommendations.
4. Silicone mat accessories should be placed in alignment with perforations in cases or trays so that sterilant pathways are not obstructed.
5. Only use internal baskets, trays, mats, or accessories that have been designed and tested for use in ETO processes. Tecomet trays and mats are made from metal, Radel® – R, or silicone and meet this requirement.

**The following table provides guidance on maximum loads. Do not overload delivery systems (cases and trays), and always follow AAMI, OSHA, and hospital standards for maximum loading. Follow instrument manufacturer’s instructions for use.**

<b>Tray Size</b>	<b>Maximum Total Weight</b>
Small (~4 x 8 x 1 in.)	1 lb. (0.45 kg)
Tray (~8 x 17 x 1 in.)	2 lb. (0.91 kg)
½ Size Case (~7 x 7 x 4 in.)	8 lb. (3.63 kg)
¾ Size Case (~7 x 12 x 4 in.)	13 lb. (5.90 kg)
Full Size Case (~9 x 19 x 4 in.)	22 lb. (9.98 kg)
Full Size Deep Vault (~9 x 19 x 6 in.)	27 lb. (12.25 kg)

## LABELING SYMBOLS

Caution

Non-Sterile

Sterile

U.S. Federal law restricts this device to sale by or on the order of a physician

CE Mark<sup>1</sup>

CE Mark with Notified Body #<sup>1</sup>

Authorized Representative in the European Community

Manufacturer

Date of Manufacture

Lot Number

Catalog Number

Consult Instructions for Use

Single Use; Do Not Reuse

Medical Device

Quantity

Country of Manufacture

<sup>1</sup>Refer to the labelling for CE information

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Symmetry Medical Inc.
Symmetry Medical Polyvac S.A.S