IMPORTANT INFORMATION - PLEASE READ BEFORE USE
INFORMATIONS IMPORTANTES - À LIRE AVANT UTILISATION
INFORMAZIONI IMPORTANTI - LEGGERE PRIMA DELL'USO
WICHTIGE INFORMATIONEN - VOR GEBRAUCH BITTE LESEN
INFORMACIÓN IMPORTANTE - LÉALA ANTES DE UTILIZAR ESTE PRODUCTO
重要な情報 - 使用前にお読みください
ВАЖНА ИНФОРМАЦИЯ - МОЛЯ, ПРОЧЕТЕТЕ ПРЕДИ УПОТРЕБА
VAŽNE INFORMACIJE - PROČITATI PRIJE UPORABE
DŮLEŽITÉ INFORMACE - PŘEČTĚTE SI PŘED POUŽITÍM
VIGTIGE OPLYSNINGER - SKAL LÆSES INDEN BRUG
BELANGRIJKE INFORMATIE - VOOR GEBRUIK DOORLEZEN
OLULINE TEAVE - PALUN LUGEGE ENNE KASUTAMIST
TÄRKEITÄ TIETOJA - LUETTAVA ENNEN KÄYTTÖÄ
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ - ΠΑΡΑΚΑΛΟΥΜΕ ΔΙΑΒΑΣΤΕ ΠΡΙΝ ΑΠΟ ΤΗ ΧΡΗΣΗ
FONTOS TUDNIVALÓK - HASZNÁLAT ELŐTT OLVASSA EL
SVARĪGA INFORMĀCIJA - LŪDZU, IZLASIET PIRMS LIETOŠANAS
SVARBI INFORMACIJA - PERSKAITYKITE PRIEŠ NAUDODAMI
VIKTIG INFORMASJON - LES FØR BRUK
WAŻNE INFORMACJE - NALEŻY PRZECZYTAĆ PRZED UŻYCIEM
INFORMAÇÕES IMPORTANTES - LER ANTES DE UTILIZAR
INFORMAȚII IMPORTANTE - VĂ RUGĂM SĂ LE CITIȚI ÎNAINTEA UTILIZĂRII
DŮLEŽITÉ INFORMÁCIE - PREČÍTAJTE SI ICH PRED POUŽITÍM
POMEMBNE INFORMACIJE - PREBERITE PRED UPORABO
VIKTIG INFORMATION - LÄS FÖRE ANVÄNDNING
ÖNEMLI BILGI - LÜTFEN KULLANMADAN ÖNCE OKUYUN

Contents (English - EN).....	1
Contenu (Français - FR).....	12
Contenuto (Italiano - IT).....	24
Inhalt (Deutsch - DE).....	36
Índice (Español - ES).....	49
内容 (日本語 - JA).....	61
Съдържание (Български - BG).....	72
Sadržaj (Hrvatski - HR).....	85
Obsah (Česky - CS).....	96
Indhold (Dansk - DA).....	107
Inhoud (Nederlands - NL).....	118
Sisukord (Eesti - ET).....	130
Sisältö (Suomi - FI).....	141
Περιεχόμενα (Ελληνικά - EL).....	153
Tartalom (Magyar - HU).....	166
Saturs (Latviski - LV).....	178
Turinys (Lietuvių - LT).....	189
Innhold (Norsk - NO).....	200
Spis treści (Polski - PL).....	211
Índice (Português - PT).....	223
Cuprins (Română - RO).....	235
Obsah (Slovenčina - SK).....	247
Vsebina (Slovenščina - SL).....	258
Innehåll (Svenska - SV).....	269
İçindekiler (Türkçe - TR).....	280

Straight Acetabular Reamer Drivers

Instructions for Use and Reprocessing

These disassembly and cleaning instructions are in accordance with ISO 17664 and AAMI ST81. They apply to:

- Reusable Straight Acetabular Reamer Drivers (provided as non-sterile)  that are intended for usage and reprocessing in a health care facility setting. All instruments and accessories may be safely and effectively reprocessed using the manual or combination manual/automated cleaning instructions and sterilization parameters provided in this document **UNLESS otherwise noted in instructions accompanying a specific instrument.**

INSTRUCTIONS FOR USE

Materials & restricted substances

For indication that the device contains a restricted substance or material of animal origin see product label.

Description

The Straight Acetabular Reamer Drivers are medical instruments they are used as part of an instrument set during a total hip replacement. The instruments are reusable and can be re-sterilized using methods defined below.

Intended Use

The Straight Acetabular Reamer Drivers are medical instruments used as part of an instrument set during a total hip replacement. The device is used to prepare the acetabulum during a total hip replacement connecting to a surgical power tool with a Tecomet Bridgeback Acetabular Reamer.

Intended Patient Population

The Straight Acetabular Reamer Drivers are prescriptive; therefore, a knowledgeable orthopedic surgeon may utilize the device on any patient they deem applicable. The device is to be used on patients undergoing a procedure that requires a total hip replacement.

Indications for Use

The Straight Acetabular Reamer Drivers are intended to prepare the acetabulum for the insertion of hip cup implants, in accordance with the implant's surgical technique and the implants indications for use and contraindications.

Contraindications

These orthopedic instruments are prescription use. The devices are only to be used by qualified health care personnel. These instruments are not for use in other surgical approaches or on other anatomical locations.

Intended User

The Straight Acetabular Reamer Drivers are prescriptive and therefore to be used by qualified orthopedic surgeons trained in the respective surgical techniques.

Accessories and/or other Device(s) intended to be used in combination with the product

- The Straight Acetabular Reamer Drivers are intended to be used with a surgical power tool via an industry standard connection.²
- The Straight Acetabular Reamer Drivers are intended to be used with Tecomet Bridgeback Acetabular Reamers.²

²For the combinations mentioned above, ensure firm connection prior to use.

Expected Clinical Benefits

When used as intended Straight Acetabular Reamer Drivers aids in the implantation of a total hip replacement cup. The implants clinical benefit is contained in the implants IFU.

Adverse Events & Complications

All surgical operations carry risks. The following are frequently encountered adverse events and complications related to having a surgical procedure in general:

- Delay to surgery caused by missing, damaged or worn instruments.
- Tissue injury and additional bone removal due to blunt, damaged, or incorrectly positioned instruments.
- Infection and toxicity due to improper processing.

Adverse events to user:

- Cuts, abrasions, contusions, or other tissue injury caused by burs, sharp edges, impaction, vibration or jamming of instruments.

Adverse Events & Complications - Reporting of Serious Incidents

Serious Incident Reporting (EU)

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, users, or other person's state of health,
- A serious public health threat.

Where further information is desired, please contact your local Tecomet sales representative. For instruments produced by another legal manufacturer, reference the manufacturer's instructions for use.

Performance and Characteristics

The Straight Acetabular Reamer Drivers are medical instruments used as part of an instrument set during a total hip replacement. The device is used to prepare the acetabulum during a total hip replacement connecting to a surgical power tool with a Tecomet Bridgeback Acetabular Reamer.

See below for assembly and disassembly instructions.

WARNINGS

-  U.S. Federal law restricts this device to sale by or on the order of a physician.
- Reusable instruments and accessories that are provided non-sterile  must be cleaned and sterilized according to these instructions prior to each use.
- The tissue protector sleeve must be disassembled for sterilization and assembled for use.
- Personal Protective Equipment (PPE) should be worn when handling, working with or disposal of contaminated or potentially contaminated instruments and accessories.
- If present, safety caps and other protective packaging material must be removed from the instruments prior to the first cleaning and sterilization.
- Caution should be exercised while handling, cleaning, wiping or disposal of instruments and accessories with sharp cutting edges, tips, and teeth.
- Ethylene oxide (EO), gas plasma, and dry heat sterilization methods are **not recommended** for sterilization of reusable instruments. Steam (moist heat) is the recommended method.

- Saline and cleaning/disinfection agents containing aldehyde, chloride, active chlorine, bromine, bromide, iodine, or iodide are corrosive and **should not** be used.
- **Do not allow biologic soil to dry on contaminated devices.** All subsequent cleaning and sterilization steps are facilitated by not allowing blood, body fluids and tissue debris to dry on used instruments.
- Automated cleaning using a washer/disinfector alone **may not** be effective for instruments with lumens, blind holes, cannulas, mated surfaces and other complex features. A thorough manual cleaning of such device features is recommended before any automated cleaning process.
- Metal brushes and scouring pads must not be used during manual cleaning. These materials will damage the surface and finish of the instruments. Use only soft bristle nylon brushes with different shapes, lengths and sizes to aid with manual cleaning.
- When processing instruments do not place heavy devices on top of delicate instruments.
- **Use of hard water should be avoided.** Softened tap water may be used for most rinsing however purified water should be used for final rinsing to prevent mineral deposits.
- Do not process instruments with polymer components at temperatures equal to or greater than 140°C/285°F because severe surface damage to the polymer will occur.
- Oils or silicone lubricants **should not** be used on surgical instruments.
- As with any surgical instrument, careful attention should be exercised to ensure that excessive force is not placed on the instrument during use. Excessive force can result in instrument failure.
- Instrument utilization is determined by the user's experience and training in surgical procedures. Do not use this instrument for any purpose outside the intended use of the device, as it may seriously affect the safety and function of the product.

Disposal

- At the end of the device's life safely dispose of the device in accordance with local procedures and guidelines.
- Any device that has been contaminated with potentially infectious substances of human origin (such as bodily fluids) should be handled according to hospital protocol for infectious medical waste. Any device that contains sharp edges should be discarded according to hospital protocol in the appropriate sharps container.

Device Life

- The Straight Acetabular Reamer Drivers is a reusable instrument. Life expectancy is dependent on the frequency of use and the care and maintenance the instruments receive. However, even with proper handling as well as correct care and maintenance, reusable instruments should not be expected to last indefinitely.
- Instruments should be inspected for damage and wear before each use. Instruments that show signs of damage or excessive wear should not be used.

Limitations on Reprocessing

- Repeated processing according to these instructions has minimal effect upon metal reusable instruments and accessories unless otherwise noted. End of life for stainless steel or other metal surgical instruments is generally determined by wear and damage incurred during the intended surgical use.
- Instruments comprised of polymers or incorporating polymer components can be sterilized using steam however they are not as durable as their metal counterparts. If polymer surfaces show signs of excessive surface damage (e.g. crazing, cracks or delamination), distortion, or visible warping they should be replaced. Contact your Tecomet representative for your replacement needs.

- Non-foaming, neutral pH enzymatic and cleaning agents are recommended for processing reusable instruments and accessories.
- Alkaline agents with a pH of 12 or less may be used to clean stainless steel and polymer instruments in countries where required by law or local ordinance or where prion diseases such as Transmissible Spongiform Encephalopathy (TSE) and Creutzfeld-Jakob Disease (CJD) are a concern. **It is critical that alkaline cleaning agents are completely and thoroughly neutralized and rinsed from the devices or degradation may occur that limits the device life.**

In countries where reprocessing requirements are more stringent than those provided in this document it is the responsibility of the user/processor to comply with those prevailing laws and ordinances.

These reprocessing instructions have been validated as being capable of preparing reusable instruments and accessories for surgical use.

 **WARNING** It is the responsibility of the user/hospital/health care provider to ensure that reprocessing is performed using the appropriate equipment, materials, and that personnel have been adequately trained in order to achieve the desired result; this normally requires that equipment and processes are validated and routinely monitored.

 **WARNING** Any deviation by the user/hospital/health care provider from these instructions should be evaluated for effectiveness to avoid potential adverse consequences.

REPROCESSING INSTRUCTIONS

Point of Use

- Remove excess biologic soil from the instruments with a disposable wipe. Place devices in a container of distilled water or cover with damp towels.

Note: Soaking in proteolytic enzyme solution prepared according to the manufacturer will facilitate cleaning especially in instruments with complex features such as lumens, mating surfaces, blind holes and cannulas.

- If instruments cannot be soaked or maintained damp, then they should be cleaned as soon as possible after use to minimize the potential for drying prior to cleaning.

Containment and Transportation

- Used instruments must be transported to the decontamination area for reprocessing in closed or covered containers to prevent unnecessary contamination risk.

Preparation for Cleaning

- Instruments designed to come apart must be **disassembled prior to cleaning and sterilization**. Disassembly of the device is described below.

Note: All recommended disassembly will be possible by hand. Never use tools to disassemble instruments beyond what is recommended.

- All cleaning solutions should be prepared at the dilution and temperature recommended by the manufacturer. Softened tap water may be used to prepare cleaning solutions.

Note: Fresh cleaning solutions should be prepared when existing solutions become grossly contaminated (turbid).

Manual Cleaning Steps

- **Step 1:** Prepare a proteolytic enzyme solution according to the manufacturer's instructions.

- **Step 2:** Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes, and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces. For device specific manipulations and maneuvering during cleaning see below.
- **Step 3:** Soak instruments for a minimum of 10 minutes. While soaking, scrub surfaces using a soft nylon-bristled brush until all visible soil has been removed. Complete Actuate moveable mechanisms. Particular attention should be given to crevices, hinged joints, box locks, instrument teeth, rough surfaces and areas with moving components or springs. Lumens, blind holes, and cannulas should be cleaned using a snug fitting round nylon bristle brush. Insert the snug fitting round brush into the lumen, blind hole, or cannula with a twisting motion while pushing in and out multiple times.

Note: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.

- **Step 4:** Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas and other difficult to access areas.
- **Step 5:** Prepare an ultrasonic cleaning bath with detergent and de-gas according to the manufacturer's recommendations. Completely submerge instruments in the cleaning solution and gently shake them to remove any trapped bubbles. Lumens, blind holes and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces. Sonically clean the instruments at the time, temperature and frequency recommended by the equipment manufacturer and optimal for the detergent used. A minimum of ten (10) minutes is recommended.

Notes:

- **For device specific manipulations and maneuvering during cleaning see below.**
- **Separate stainless steel instruments from other metal instruments during ultrasonic cleaning to avoid electrolysis.**
- **Fully open hinged instruments.**
- **Use wire mesh baskets or trays designed for ultrasonic cleaners.**
- **Regular monitoring of sonic cleaning performance by means of an ultrasonic activity detector, aluminum foil test, TOSI™ or SonoCheck™ is recommended.**
- **Step 6:** Remove the instruments from the ultrasonic bath and rinse in purified water for a minimum of one (1) minute or until there is no sign of residue detergent or biologic soil. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas and other difficult to access areas.
- **Step 7:** Dry instruments with a clean, absorbent non-shedding lint free cloth. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas, and difficult to access areas.

Combination Manual/Automated Cleaning Steps

- **Step 1:** Prepare a proteolytic enzyme solution according to the manufacturer's instructions.
- **Step 2:** Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes and cannulations

should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces. For manipulations and maneuvering during cleaning see below.

- **Step 3:** Soak instruments for a minimum of 10 minutes. Scrub surfaces using a soft nylon-bristled brush until all visible soil has been removed. Actuate moveable mechanisms. Particular attention should be given to crevices, hinged joints, box locks, instrument teeth, roughened surfaces, and areas with moving components or springs. Lumens, blind holes and cannulas should be cleaned using a snug fitting round nylon bristle brush. Insert the snug fitting round nylon bristle brush into the lumen, blind hole or cannula with a twisting motion while pushing in and out multiple times.

Notes:

- **For device specific manipulations and maneuvering during cleaning see below.**
- **All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.**
- **Step 4:** Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas and other difficult to access areas.
- **Step 5:** Place instruments in a suitable validated washer/ disinfecter. Follow the washer/disinfecter manufacturer's instructions for loading the instruments for maximum cleaning exposure, e.g. open all instruments, place concave instruments on their side or upside down, use baskets and trays designed for washers, place heavier instruments on the bottom of trays and baskets. If the washer/disinfecter is equipped with special racks (e.g. for cannulated instruments) use them according to the manufacturer's instructions.
- **Step 6:** Process instruments using a standard washer/ disinfecter instrument cycle according to the manufacturer's instructions. The following minimum wash cycle parameters are recommended:

Cycle	Description
1	Pre-wash ● Cold Softened Tap Water ● 2 minutes
2	Enzyme Spray & Soak ● Hot Softened Tap Water ● 1 minute
3	Rinse ● Cold Softened Tap Water
4	Detergent Wash ● Hot Tap Water (64-66°C/146-150°F) ● 2 minutes
5	Rinse ● Hot Purified Water (64-66°C/146-150°F) ● 1 minute
6	Hot Air Dry (116°C/240°F) ● 7 - 30 minutes

Notes:

- **The washer/disinfecter manufacturer's instructions should be followed.**
- **A washer/disinfecter with demonstrated efficacy (e.g. FDA clearance or validated to ISO 15883) should be used.**
- **Dry time is shown as a range because it is dependent upon the load size placed into the washer/disinfecter.**
- **Many manufacturers pre-program their washer/ disinfecters with standard cycles and they may include a thermal low-level disinfection rinse after the detergent wash. The thermal disinfection cycle should be performed to achieve a minimum value A0 = 600 (e.g. 90°C/194°F for 1 minute according to ISO 15883-1) and is compatible with the instruments.**
- **If a lubrication cycle is available that applies to water-soluble lubricant such as Preserve®, Instrument Milk, or equivalent material intended for medical device application**

it is acceptable to use on instruments unless otherwise indicated.

Disinfection

- Instruments must be terminally sterilized prior to use. See sterilization instructions below.
- Low level disinfection may be used as part of a washer/disinfectant cycle, but the devices must also be sterilized before use.

Drying

- Dry instruments with a clean, absorbent non-shedding lint free cloth. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas and difficult to access areas.

Inspection & Testing

- After cleaning, all devices should be thoroughly inspected for residue biologic soil or detergent. If contamination is still present repeat the cleaning process.
- Visually inspect each device for completeness, damage, and excessive wear. If damage or wear is observed that might compromise the function of the device, do not process them further and contact your Tecomet representative for a replacement.
- When inspecting devices look for the following:
 - Cutting edges should be free of nicks and have a continuous edge.
 - Jaws and teeth should align properly.
 - Movable parts should operate smoothly throughout the intended range of motion.
 - Locking mechanisms should fasten securely and close easily.
 - Long thin instruments should be free of bending or distortion.
 - Where instruments form part of a larger assembly, check that all components are available and assemble readily.
 - Polymer surfaces should not show signs of excessive surface damage (e.g. crazing, cracks or delamination), distortion, or visible warping. If the instrument is damaged, it should be replaced.

Lubrication

- After cleaning and before sterilization, instruments with moving parts (e.g. hinges, box-locks, sliding or rotating parts) should be lubricated with a water-soluble lubricant such as Preserve®, Instrument Milk, or equivalent material intended for medical device application. Always follow the lubricant manufacturer's instructions for dilution, shelf life and application method.

Packaging for Sterilization

- Instruments and container must be properly cleaned before sterilization.
- Place the instruments in their respective position within the delivery system according to the markings/labeling in the container.
- Once the container is loaded, put the lid on and secure all latches or locks.
- Single devices may be packaged in an approved (e.g. FDA cleared or ISO 11607 compliant) medical grade sterilization pouch or wrap. Care should be used when packaging so that the pouch or wrap is not torn. Devices should be wrapped using the double wrap or equivalent method (ref: AAMI ST79, AORN Guidelines).
- Reusable wraps are not recommended.
 - The case or tray must be wrapped in an approved (e.g. FDA cleared or ISO 11607 compliant) medical grade sterilization wrap by following the double wrap method or equivalent (ref: AAMI ST79, AORN Guidelines).
 - Follow the case/tray manufacturer's recommendations for loading and weight. Total weight of the wrapped case or tray should not exceed 11.4kg/25lbs.

Sterilization

- Moist heat/steam sterilization is the recommended method for the instruments.
- Use of an approved chemical indicator (class 5) or chemical emulator (class 6) within each sterilization load is recommended.
- Always consult and follow the sterilizer manufacturer instructions for load configuration and equipment operation. Sterilizing equipment should have demonstrated efficacy (e.g. FDA clearance, EN 13060 or EN 285 compliance). Additionally, the manufacturer's recommendations for installation, validation, and maintenance should be followed.
- Validated exposure times and temperatures to achieve a 10⁻⁶ sterility assurance level (SAL) are listed in the following table.
- Local or national specifications should be followed where steam sterilization requirements are **more conservative** than those listed in the table below.

Cycle Type	Temperature	Exposure Time	Dry Time
United States Recommended Parameters			
Pre-vacuum / Vacuum Pulse	132°C / 270°F	4 minutes	40 minutes
Cycle Type	Temperature	Exposure Time	Dry Time
European Recommended Parameters			
Pre-vacuum/ Vacuum Pulse	134°C / 273°F	3 minutes	40 minutes

Drying & Cooling

- The recommended drying time for single wrapped instruments is 30 minutes unless otherwise noted in device specification instructions.
- A 15 minute open door time is recommended after drying.
- A 30 minute minimum cooling time is recommended after drying but longer times may be necessary because of load configuration, ambient temperature and humidity, device design, and packaging used.

Note: Disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is a concern about TSE/CJD contamination are: 134°C/273°F for 18 minutes. These devices are compatible with these parameters.

Storage

- Sterile packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes.

Note: Inspect every package before use to ensure that the sterile barrier (e.g. wrap, pouch or filter) is not torn, perforated, showing signs of moisture, or appearing to be tampered with. If any of those conditions are present then the contents are considered non-sterile and should be re-processed through cleaning, packaging and sterilization.


Families of Straight Acetabular Reamer Drivers

The Straight Acetabular Reamer Drivers come in two styles: Bridgeback and Dual Connection. Both the Bridgeback and the Dual Connection Acetabular Reamer Drivers have a quick connection that releases the Acetabular Reamer through use of a Release Collar. Furthermore, they both have a Tissue Protector, a spring that loads the quick connection, and a spring-loaded Release Button as seen in Figure 1.

Two styles of Tissue Protectors are offered: a fully retained version and a fully removable as shown in Figure 1. With the fully removable Tissue Protector, it is important that the spring and Collar are retained on the shaft of the Acetabular Reamer Driver. The Acetabular Reamer Drivers come in various lengths, so it is

important to assemble the correct length Tissue Protector with the appropriate shaft.

Figure 1. Components of an Acetabular Reamer Driver (Bridgeback is shown)

Dual Connection Straight Acetabular Reamer Driver Disassembly

1. Remove the Acetabular Reamer from the Acetabular Reamer Driver by pulling the Release Collar proximally as shown in Figure 2.
2. Remove the Tissue Protector proximally by depressing the release button and sliding the Tissue Protector proximally as shown in Figures 3 and 4.

Figure 2. Disconnection of Acetabular Reamer from Dual Connection Driver


Figure 3. Proximal Movement of Tissue Protector on Dual Connection Driver



Figure 4. Remove Tissue Protector for Cleaning and Sterilization on Dual Connection Driver

Standard and Extended Length Acetabular Reamer Driver (Bridgeback) Disassembly

1. Remove the Acetabular Reamer from the Acetabular Reamer Driver by pulling the Release Collar proximally and rotating the Acetabular Reamer as shown in Figure 5.
2. Remove the Tissue Protector by depressing the Release Button and sliding the Tissue Protector proximally as shown in Figure 6. The spring and Release Collar should remain on the Acetabular Reamer Driver shaft.

Figure 5. Disconnection of the Acetabular Reamer from the Bridgeback Driver



Figure 6. Disassembly of Tissue Protector on Bridgeback Reamer Driver

Manipulations and Maneuvering During Cleaning



1. While immersed, actuate and maneuver all moving parts of the instruments to ensure complete exposure of the cleaning solution to all surfaces.
2. Depress Release Button multiple times while cleaning.
3. Move spring and Release Collar proximally and scrub the distal tip of the Acetabular Reamer Driver.
 - a. Ensure that the spring is maintained on the shaft.
4. For removable Tissue Protectors, scrub the inside of the Tissue Protector with a snug fitting pipe cleaner. Pass the pipe cleaner through the inside of the tube three (3) times.
5. Flush the inside of the Tissue Protector with a syringe.
6. Rinse all parts thoroughly with distilled or sterile water until all traces of cleaning solution are removed. Maneuver all moving parts while rinsing.
7. Lubricate the moving parts of the Acetabular Reamer Driver with water soluble lubricant that is intended for surgical instruments.

Assembly Prior to Use for Drivers with Removable Tissue Protectors


Prior to use, ensure that the device is appropriately assembled.

1. Ensure that the correct length Tissue Protector is utilized with the appropriate Acetabular Reamer Driver.
 - For the Bridgeback design, check that the etch description on the Tissue Protector matches the etched description on the Acetabular Reamer Driver.
2. The spring must be present and functional.
3. The stepped end of the Tissue Protector must go on the Acetabular Reamer Driver first (see Figure 7).
4. Depress the Release Button and slide the Tissue Protector over it.
5. Assess fit of the Tissue Protector.
 - There should not be a lot of excess space between the Tissue Protector and the Release Button.
 - If you believe there is excess space on the Bridgeback design, remove the Tissue Protector and ensure that the etched length (Standard or Extended) on the Acetabular Reamer Driver match the length on the Tissue Protector.

Figure 7. Stepped End of the Tissue Protectors

SYMBOLS USED ON LABELING

-  Caution.
-  Lot Number
-  Non-Sterile
-  Catalog Number
-  U.S Federal law restricts this device to sale by or on the order of a physician
-  Consult Instructions for Use
-  CE Mark with Notified Body.¹
-  Medical Device
-  Authorized Representative in European Community
-  Packaging Unit
-  Swiss Authorized Representative
-  Country of Manufacture
-  Manufacturer
-  Distributor
-  Date of Manufacture
-  Use By
-  Importer
-  Unique Device Identifier

¹Refer to the labeling for CE information

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Symmetry Medical Manufacturing, Inc.
Symmetry Medical Polyvac S.A.S
DePuy Orthopaedics, Inc