IMPORTANT INFORMATION - PLEASE READ BEFORE USE
INFORMATIONS IMPORTANTES - À LIRE AVANT UTILISATION
INFORMAZIONI IMPORTANTI - LEGGERE PRIMA DELL'USO
WICHTIGE INFORMATIONEN - VOR GEBRAUCH BITTE LESEN
INFORMACIÓN IMPORTANTE - LÉALA ANTES DE UTILIZAR ESTOS PRODUCTOS
重要な情報 - 使用前にお読みください
ВАЖНА ИНФОРМАЦИЯ - МОЛЯ, ПРОЧЕТЕТЕ ПРЕДИ УПОТРЕБА
VAŽNE INFORMACIJE - PROČITATI PRIJE UPORABE
DŮLEŽITÉ INFORMACE - PŘEČTĚTE SI PŘED POUŽITÍM
VIGTIGE OPLYSNINGER - LÆS VENLIGST INDEN BRUG
BELANGRIJKE INFORMATIE - VOOR GEBRUIK DOORLEZEN
OLULINE TEAVE - PALUN LUGEGE ENNE KASUTAMIST
TÄRKEITÄ TIETOJA - LUETTAVA ENNEN KÄYTTÖÄ
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ -
ΠΑΡΑΚΑΛΟΥΜΕ ΔΙΑΒΑΣΤΕ ΠΡΙΝ ΑΠΟ ΤΗ ΧΡΗΣΗ
FONTOS TUDNIVALÓK - HASZNÁLAT ELŐTT OLVASSA EL
SVARĪGA INFORMĀCIJA - LŪDZU, IZLASIET PIRMS LIETOŠANAS
SVARBI INFORMACIJA - PRIEŠ NAUDOJANTIS BŪTINA PERSKAITYTI
VIKTIG INFORMASJON - LES FØR BRUK
WAŽNE INFORMACJE - NALEŻY PRZECZYTAĆ PRZED UŻYCIEM
INFORMAÇÃOES IMPORTANTES - LER ANTES DE UTILIZAR
INFORMAȚII IMPORTANTE - VĂ RUGĂM SĂ LE CITIȚI ÎNAINTEA UTILIZĂRII
DŔLEŽITÉ INFORMÁCIE - PREČÍTAJTE SI PRED POUŽITÍM
POMEMBNE INFORMACIJE - PREBERITE PRED UPORABO
VIKTIG INFORMATION - LÄS FÖRE ANVÄNDNING
ÖNEMLİ BİLGİ - LÜTFEN KULLANMADAN ÖNCE OKUYUN

Contents (English - EN)	1
Contenu (Français - FR)	6
Contenuto (Italiano - IT)	11
Inhalt (Deutsch - DE)	16
Índice (Español - ES)	21
内容 (日本語 - JA)	26
Съдържание (Български - BG)	30
Sadržaj (Hrvatski - HR)	35
Obsah (Česky - CS)	40
Indhold (Dansk - DA)	45
Inhoud (Nederlands - NL)	50
Sisukord (Eesti - ET)	55
Sisältö (Suomi - FI)	60
Περιεχόμενα (Ελληνικά - EL)	65
Tartalom (Magyar - HU)	71
Saturs (Latviski - LV)	76
Turinys (Lietuvių - LT)	81
Innhold (Norsk - NO)	86
Spis treści (Polski - PL)	91
Índice (Português - PT)	96
Cuprins (Română - RO)	101
Obsah (Slovenčina - SK)	106
Vsebina (Slovenščina - SL)	111
Innehåll (Svenska - SV)	116
İçindekiler (Türkçe - TR)	121

Instructions For Use and Reprocessing for Tecomet Hybrid Cases and Trays

Product Description/Intended Use

The Tecomet Hybrid (Metal Lid/Base, Radel Inserts) Cases and Trays are intended for the protection, organization, sterilization, and delivery to the instruments to the surgical field. The Case and Trays are not intended to maintain sterility by themselves. They are designed to facilitate the sterilization process when used with an FDA approved sterilization wrap. The wrapping material is designed to allow air removal, steam penetration/evacuation (drying) and to maintain the sterility of the internal components.

Indications For Use

The Tecomet Hybrid Cases and Trays are indicated for use to protect, organize, transport, and facilitate the sterilization process of the intended surgical instruments by allowing steam penetration and air removal. When used in conjunction with an approved sterilization wrap and following the specified IFU, the sterility of the enclosed medical device is maintained provided the sterilized instruments are properly housed and sterility maintenance is performed per ANSI ST79.

Contraindications

The Tecomet Hybrid Cases and Trays have **NOT** been validated for use with flexible endoscopes, or devices with lumens or working channels longer than 4 inches (3mm ID). Always refer to instrument manufacturer instructions.

Intended Patient Population

None. The Tecomet Hybrid Cases and Trays are intended to protect, organize, transport, and facilitate the sterilization of the intended surgical instruments.

Intended User

The Tecomet Hybrid Cases and Trays are intended to be used by qualified healthcare professionals, such as surgeons, physicians, operating room staff, and individuals involved in the preparation of the intended surgical instruments.

Expected Clinical Benefits

The Tecomet Hybrid Cases and Trays protect, organize, transport, and facilitate the sterilization of the intended surgical instruments.

Performance Characteristics

The Tecomet Hybrid Cases and Trays protect, organize, transport, and facilitate the sterilization of the intended surgical instruments. The performance and safety of the cases are established, representing the current state of the art when used as intended.

Adverse Events & Complications

All surgical operations carry risk. The following are frequently encountered adverse events and complications related to having a surgical procedure in general:

- Delay to surgery caused by missing, damaged or worn instruments.
- Tissue injury and additional bone removal due to blunt, damaged or incorrectly positioned instruments.
- Infection and toxicity due to improper processing.

Adverse events to user:

- Cuts, abrasions, contusions or other tissue injury caused by burrs, sharp edges, impaction, vibration or jamming of instruments.

Adverse Events & Complications – Reporting of Serious Incidents

Serious Incident Reporting (EU)

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

Materials and Restricted Substance

For indication that the device contains a restricted substance or material of animal origin see product label.

Combination of Medical Devices

The Tecomet Hybrid Cases and Trays are designed to protect, organize, transport, and facilitate the sterilization of the intended surgical instruments. The Hybrid Cases and Trays are not to be used in combination with instruments provided by other manufacturers. Tecomet assumes no liability in such instances.

Disposal

- At the end of the device's life safely dispose of the device in accordance with local procedures and guidelines.
- Any device that has been contaminated with potentially infectious substances of human origin (such as bodily fluids) should be handled according to hospital protocol for infectious medical waste. Any disposable or defective device that contains sharp edges should be discarded according to hospital protocol in the appropriate sharps container.

WARNING

-  U.S Federal law restricts this device to sale by or on the order of a physician.
- The cases are provided NON-STERILE  and must be cleaned and sterilized according to these instructions prior to each use.
- Personal Protective Equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated cases.
- Saline and cleaning/disinfection agents containing aldehyde, chloride, active chlorine, bromine, bromide, iodine or iodide are corrosive and **should not** be used.
- **Do not allow biologic soil to dry on contaminated devices.** All subsequent cleaning and sterilization steps are facilitated by not allowing blood, body fluids and tissue debris to dry on used cases.
- Metal brushes and scouring pads **MUST NOT** be used during manual cleaning. These materials will damage the surface and finish of the cases. Use only soft bristle nylon brushes with different shapes, lengths and sizes to aid with manual cleaning.
- **Use of hard water should be avoided.** Softened tap water may be used for most rinsing however purified water should be used for final rinsing to prevent mineral deposits.

Cleaning

- Both physical and chemical (detergent) processes may be necessary to clean soiled items.
- Chemical (detergent) cleaners alone cannot remove all soil and debris; therefore, a careful manual cleaning of each item with a soft sponge or cloth is essential for maximum decontamination.
- For difficult areas, a clean soft bristled brush is recommended.
- Once the items have been cleaned, they should be thoroughly rinsed with clean water to remove any detergent or chemical residue before sterilization.

- **Tecomet** recommends the use of a mild enzymatic detergent with a near neutral pH.
- Do not use solvents, abrasive cleaners metal brushes, or abrasive pads.
- Case and trays may be placed in mechanical cleaning equipment.

Packaging for Sterilization

- Reference Case markings for placing the cleaned instruments to their respective locations.
- **Always** inspect instruments and case for cleanliness and damage before use. Make sure all latches and handles are secure and in working order. Do not overload cases and trays. Balance contents uniformly within the container and arrange to allow steam to come into contact with all objects in the container.
- Tecomet recommends that the cases and trays be wrapped according to the sterilization wrap manufacturer’s instructions prior to sterilization to maintain sterility of internal components/ items and for proper aseptic presentation to the surgical field.
- Testing was performed with the Acetabular Reamer Case and the instruments included in the kit. Always follow instrument manufacturer instructions if their sterilization or drying recommendations exceed these guidelines.

Sterilization

- Moist heat/steam sterilization is the recommended method for the device.
- Use of an approved chemical indicator (class 5) or chemical emulator (class 6) within each sterilization load is recommended.
- Always consult and follow the sterilizer manufacturer instructions for load configuration and equipment operation. Sterilizing equipment should have demonstrated efficacy (e.g., FDA clearance, EN 13060 or EN 285 compliance,). Additionally, the manufacturer’s recommendations for installation, validation, and maintenance should be followed.
- Validated exposure times and temperatures to achieve a 10⁻⁶ sterility assurance level (SAL) are listed in the following table.
- Wrapped cases, trays, and instruments should be sterilized according to the table below.

Cycle Type	Temperature	Exposure Time	Dry Time	Cool-Down Time*
United States Recommended Parameters				
Pre-vacuum / Vacuum Pulse	132°C / 270°F	4 minutes	40 minutes	30 minutes
Cycle Type	Temperature	Exposure Time	Dry Time	Cool-Down Time*
European Recommended Parameters				
Pre-vacuum/ Vacuum Pulse	134°C / 273°F	3 minutes	40 minutes	30 minutes

*Allow the wrapped device to cool outside of the autoclave for 30 minutes prior to use.

Filtered Sterilization Container Systems

The Tecomet Hybrid Cases and Trays have not been validated for use in filtered sterilization containers and Tecomet does not recommend use of such a system. If filtered sterilization container systems are used, the user is responsible for following the recommendations of the manufacturer for the proper placement and use of cases and tray inside the container.

Ethylene Oxide Sterilization (EtO)

- The Tecomet Hybrid Cases and Trays have **NOT** been validated for EtO sterilization of devices.

Notice

1. **DO NOT** load cases into sterilizer on sides or upside down with lid side on shelf or cart. Load cases on cart or shelf so that the lid is always facing upwards. This will allow for proper drying. The hybrid cases and trays are designed to drain in this position.
2. **DO NOT STACK** cases or trays in an autoclave chamber. Stacking the cases and overloading the cases will adversely affect sterilization and drying effectiveness.
3. After the autoclave door is opened, all cases must be allowed to cool thoroughly. Place cases on rack or shelf with linen cover until cooling is complete. The potential for condensation may increase if the case is not allowed to cool properly.
4. If condensation is observed, check to ensure steps 1, 2, and 3 were followed. In addition, verify that the steam, which is being used for sterilization processing, has a quality of more than 97%¹. Also confirm that the sterilizers have been inspected for routine maintenance in accordance with manufacturer recommendations.
5. Silicone mat accessories should be placed in alignment with the perforations in the case or tray so that the sterilant pathways are not obstructed.

The following table provides guidance on maximum loads. Do not overload delivery systems (cases and trays), and always follow AAMI, OSHA and hospital standards for maximum loading. Follow instrument manufacturer's instructions for use.

Tray Size	Maximum Total Weight
Small (~ 4 x 8 x 1 in.)	1 lb. (0.45 kg)
Tray (~ 8 x 17 x 1 in.)	2 lb. (0.91 kg)
½ Size Case (~ 7 x 7 x 4 in.)	8 lb. (3.6 kg)
¾ Size Case (~ 7 x 12 x 4 in.)	13 lb. (5.9 kg)
Full Size Case (~ 9 x 19 x 4 in.)	22 lb. (10.0 kg)

REFERENCES

¹ANSI/AAIM ST79: Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities

SYMBOLS USED ON LABELING

Caution

Non-Sterile

U.S Federal law restricts this device to sale by or on the order of a physician

CE Mark

Authorized Representative in the European Community

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Symmetry Medical Manufacturing, Inc.
Symmetry Medical Polyvac S.A.S