IMPORTANT INFORMATION - PLEASE READ BEFORE USE  
INFORMATIONS IMPORTANTES - À LIRE AVANT UTILISATION  
INFORMAZIONI IMPORTANTI - LEGGERE PRIMA DELL'USO  
WICHTIGE INFORMATIONEN - VOR GEBRAUCH BITTE LESEN  
INFORMACIÓN IMPORTANTE - LÉALA ANTES DE UTILIZAR ESTE PRODUCTO  
重要な情報 - 使用前にお読みください  
ВАЖНА ИНФОРМАЦИЯ - МОЛЯ, ПРОЧЕТЕТЕ ПРЕДИ УПОТРЕБА  
VAŽNE INFORMACIJE - PROČITATI PRIJE UPORABE  
DŮLEŽITÉ INFORMACE - PŘEČTĚTE SI PŘED POUŽITÍM  
VIGTIGE OPLYSNINGER - SKAL LÆSES INDEN BRUG  
BELANGRIJKE INFORMATIE - VOOR GEBRUIK DOORLEZEN  
OLULINE TEAVE - PALUN LUGEGE ENNE KASUTAMIST  
TÄRKEITÄ TIETOJA - LUETTAVA ENNEN KÄYTTÖÄ  
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ - ΠΑΡΑΚΑΛΟΥΜΕ ΔΙΑΒΑΣΤΕ ΠΡΙΝ ΑΠΟ ΤΗ ΧΡΗΣΗ  
FONTOS TUDNIVALÓK - HASZNÁLAT ELŐTT OLVASSA EL  
SVARĪGA INFORMĀCIJA - LŪDZU, IZLASIET PIRMS LIETOŠANAS  
SVARBI INFORMACIJA - PERSKAITYKITE PRIEŠ NAUDODAMI  
VIKTIG INFORMASJON - LES FØR BRUK  
WAŻNE INFORMACJE - NALEŻY PRZECZYTAĆ PRZED UŻYCIEM  
INFORMAÇÃOES IMPORTANTES - LER ANTES DE UTILIZAR  
INFORMAȚII IMPORTANTE - VĂ RUGĂM SĂ LE CITIȚI ÎNAINTEA UTILIZĂRII  
DÔLEŽITÉ INFORMÁCIE - PREČÍTAJTE SI ICH PRED POUŽITÍM  
POMEMBNE INFORMACIJE - PREBERITE PRED UPORABO  
VIKTIG INFORMATION - LÄS FÖRE ANVÄNDNING  
ÖNEMLİ BİLGİ - LÜTFEN KULLANMADAN ÖNCE OKUYUN

Contents (English – EN).....	1
Contenu (Français – FR).....	6
Contenuto (Italiano – IT).....	12
Inhalt (Deutsch – DE).....	17
Índice (Español – ES).....	22
内容 (日本語 – JA).....	27
Съдържание (български – BG).....	32
Sadržaj (Hrvatski – HR).....	38
Obsah (Česky – CS).....	43
Indhold (Dansk – DA).....	48
Inhoud (Nederlands – NL).....	53
Sisukord (Eesti – ET).....	58
Sisältö (Suomi – FI).....	63
Περιεχόμενα (Ελληνικά – EL).....	68
Tartalom (Magyar – HU).....	74
Saturs (Latviski – LV).....	79
Turinys (Lietuvių kalba – LT).....	84
Innhold (Norsk – NO).....	89
Spis treści (Polski – PL).....	94
Índice (Português – PT).....	99
Cuprins (Română – RO).....	104
Obsah (Slovenčina – SK).....	109
Vsebina (Slovenščina – SL).....	114
Innehåll (Svenska – SV).....	119
İçindekiler (Türkçe – TR).....	124

### INSTRUCTIONS FOR USE/INTENDED USE

Tecomet® delivery systems (cases and trays) are intended for the protection, organization and the delivery to the surgical field of surgical instruments and/or other medical devices.

Tecomet delivery systems (cases and trays) are not designed to maintain sterility by themselves. They are designed to facilitate the sterilization process when used in conjunction with a wrapping material (FDA cleared sterilization wrap) or a specified filtered sterilization container system. Wrapping materials and sterilization containers are designed to allow air removal, steam penetration/evacuation (drying) and maintain the sterility of the internal components.

### Recommendations for Care, Cleaning and Sterilization of Tecomet Delivery Systems (cases and trays)

Both physical and chemical (detergent) processes may be necessary to clean soiled items. Chemical (detergent) cleaners alone cannot remove all soil and debris; therefore, a careful manual cleaning of each item with a soft sponge or cloth is essential for maximum decontamination. For difficult access areas, a clean soft bristled brush is recommended. Once the items have been cleaned, they should be thoroughly rinsed with clean water to remove any detergent or chemical residue before sterilization. **Tecomet** recommends the use of a mild enzymatic detergent with a near neutral pH. Do not use solvents, abrasive cleaners, metal brushes or abrasive pads. Cases and trays may be placed in mechanical cleaning equipment.

**Always** inspect for cleanliness or damage before use. Make sure all latches and handles are secure and in working order. Do not overload cases and trays. Balance contents uniformly within the container and arrange to allow steam to come in contact with all objects in the container. Tecomet recommends that the delivery systems (cases and trays) be processed according to the sterilization wrap or filtered container manufacturers instructions prior to sterilization to maintain sterility of internal components/ items and for proper aseptic presentation to the surgical field.

For steam autoclaves, Tecomet has verified product performance in the following cycles. Our testing was done with common surgical appliances such as reamers, drills, hammers, rasps, drivers, chisels, awls, handpieces, rongeurs, blades and bits, including short lumened devices.

Complex instruments may require disassembly and extended sterilization times. Always follow instrument manufacturer instructions if their sterilization or drying recommendations exceed these guidelines.

### **1. Prevacuum Sterilizer**

Wrapped cases, trays and instruments, or cases, trays and instruments should be exposed to 132° C to 135° C (270° F to 275° F) for at least four minutes. Dry for 20 to 40 minutes.

### **2. Gravity Displacement Sterilizer**

Wrapped cases, trays and instruments should be exposed to 132° C to 135° C (270° F to 275° F) for at least 30 minutes, or 121° C to 123° C (250° F to 254° F) for at least 55 minutes. Dry for 20 to 50 minutes.

Variables that may affect drying times include: loading density of the case/tray, instrument configuration, total contents of the sterilizer, steam quality, equipment maintenance and others.

### **3. Filtered Sterilization Container Systems**

Tecomet delivery systems (cases and trays) have been validated for placement inside Ultra Container, Genesis and Aesculap filtered sterilization container systems and exposed to 132° C to 135° C (270° F to 275° F) for at least four minutes for prevacuum sterilizers only. Tecomet delivery systems have not been validated for use in filtered sterilization containers in gravity displacement sterilizers and should not be used in such a system.

Follow the recommendations of the manufacturer of the filtered sterilization containers for the proper placement and use of cases and trays inside the container.

### **Contraindications**

Tecomet delivery systems (cases and trays) have NOT been validated for ETO sterilization of devices with lumens or working channels longer than 8-1/2 inches (1.5mm ID). Always refer to instrument manufacturer instructions.

Stacking of delivery systems (cases and trays) and overloading of the units will adversely effect sterilization and drying effectiveness.

**DO NOT STACK** cases and trays in autoclave chamber.

**Notice:**

1. **DO NOT** load cases into sterilizer on sides or upside down with lid side on the shelf or cart. Load cases on cart or shelf, so that the lid is always facing upward. This will allow for proper drying. Tecomet cases and trays are designed to drain in this position.
2. After the autoclave door is opened, all cases must be allowed to cool thoroughly. Place cases on a rack or shelf with linen cover until cooling is complete. The potential for condensation may increase if the case is not allowed to cool properly.
3. If condensation is observed, check to ensure steps 1 and 2 were followed. In addition, verify that the steam, which is being used for sterilization processing, has a quality of more than 97%. Also confirm that the sterilizers have been inspected for routine maintenance in accordance with manufacturer recommendations.
4. Silicone mat accessories should be placed in alignment with perforations in cases or trays so that sterilant pathways are not obstructed.
5. Only use internal baskets, trays, mats, or accessories that have been designed and tested for use in ETO processes. Tecomet trays and mats are made from metal, Radel<sup>®</sup>-R or silicone and meet this requirement.

## Symbols Used on Labeling

Manufacturer

Catalog Number

Caution. Consult Accompanying Documents

Consult Instructions for Use

Lot Number

Non-Sterile Instrument

U.S. Federal law restricts this device to sale by or on the order of a physician

CE Mark

Authorized Representative in the European Community

Date of Manufacture

Medical Device

Packaging Unit

Country of Manufacture

Distributor

**The following table provides guidance on maximum loads. Do not overload delivery systems (cases and trays), and always follow AAMI, OSHA and hospital standards for maximum loading. Follow instrument manufacturers instructions for use.**

<b>Tray Size</b>	<b>Maximum Total Weight</b>
Small (~ 4 x 8 x 1 in.)	1 lb. (0.45 Kg)
Tray (~ 8 x 17 x 1 in.)	2 lb. (0.91 Kg)
1/2 Size Case (~ 7 x 7 x 4 in.)	8 lb. (3.6 Kg)
3/4 Size Case (~ 7 x 12 x 4 in.)	13 lb. (5.9 Kg)
Full Size Case (~ 9 x 19 x 4 in.)	22 lb. (10 Kg)
Full Size Deep Vault ( ~ 9 x 19 x 6 in.)	27 lb. (12.2 Kg)

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All rights reserved.
Symmetry Medical Polyvac S.A.S
KARL STORZ Endoscopy-America, Inc.
KARL STORZ Endoscopy Canada, Ltd.
Printed in USA. ©2021 Tecomet, Inc.