IMPORTANT INFORMATION - PLEASE READ BEFORE USE  
INFORMATIONS IMPORTANTES - À LIRE AVANT UTILISATION  
INFORMAZIONI IMPORTANTI - LEGGERE PRIMA DELL'USO  
WICHTIGE INFORMATIONEN - VOR GEBRAUCH BITTE LESEN  
INFORMACIÓN IMPORTANTE - LÉALA ANTES DE UTILIZAR ESTE PRODUCTO  
重要な情報 - 使用前にお読みください

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# **Evolution-TL and Evolution-C Removal**

**EN**

## **Kits Instructions for Use**

These instructions are in accordance with ISO17664 and AAMI ST81. They apply to reusable surgical instruments and accessories (provided as non-sterile) supplied by Tecomet and intended for reprocessing in a health care facility setting. All instruments and accessories may be safely and effectively reprocessed using the manual or combination manual/automated cleaning instructions and sterilization parameters provided in this document **UNLESS otherwise noted in instructions accompanying a specific instrument.**

## **Description**

The reusable manual orthopedic surgical instruments contained in the Evolution Removal Kits consist of drivers for standard metric and English hex drivers, hexalobe drivers, broken and stripped screw removers, specific drivers, extensions, ratchet handle, multi-axial handle, counter torque wrench, combination wrenches, needle-nose locking pliers and cutting instruments. The instruments are provided in an instrument storage container.

## **Intended Use**

Product Family/Device Description	Intended Use	Indication For Use	Patient Group Target	Performance Characteristics of the Device	Expected Clinical Benefits
Evolution-TL Spinal Hardware Removal Kit	Remove spinal hardware from the lumbar and thoracic regions of the spine.	Any surgical scenario that requires the removal of implants in the thoracolumbar region of the spine.	Any patient that requires the removal of implanted spinal hardware.	The Evolution-TL Spinal Hardware Removal Kits are designed to aid in the removal of implanted spinal hardware.	When used as intended, the Evolution-TL Spinal Hardware Removal Kits aid in the removal of implanted spinal hardware.

<b>Product Family/Device Description</b>	<b>Intended Use</b>	<b>Indication For Use</b>	<b>Patient Group Target</b>	<b>Performance Characteristics of the Device</b>	<b>Expected Clinical Benefits</b>
Evolution-C Spinal Removal Kit	Remove spinal hardware from the cervical region of the spine.	Any surgical scenario that requires the removal of implants in the cervical region of the spine.	Any patient that requires the removal of implanted spinal hardware.	The Evolution-C Spinal Hardware Removal Kits are designed to aid in the removal of implanted spinal hardware.	When used as intended, the Evolution-C Spinal Hardware Removal Kits aid in the removal of implanted spinal hardware.
Male Hexalobe Driver	Remove implanted screws that have a female Hexalobe configuration.	Any surgical scenario that requires the removal of an implanted screw.	Any patient that requires the removal of an implanted screw.	The Male Hexalobe Drivers are designed to aid in the removal of implanted screws that have a female Hexalobe configuration.	When used as intended, the Male Hexalobe Drivers, aid in the removal of implanted screws that have a female Hexalobe configuration.

<b>Product Family/Device Description</b>	<b>Intended Use</b>	<b>Indication For Use</b>	<b>Patient Group Target</b>	<b>Performance Characteristics of the Device</b>	<b>Expected Clinical Benefits</b>
Male Hex Driver	Remove implanted screws that have a female hex configuration.	Any surgical scenario that requires the removal of an implanted screw.	Any patient that requires the removal of an implanted screw.	The Male Hex Drivers are designed to aid in the removal of implanted screws that have a female hex configuration.	When used as intended, the Male Hex Drivers, aid in the removal of implanted screws that have a female hex configuration.
Ratcheting Driver	Increase efficiency when tightening or loosening fasteners.	Any surgical scenario where a clinician desires to increase efficiency when tightening or loosening fasteners.	Any patient that requires tightening or loosening fasteners.	The Ratcheting Driver are designed to aid in tightening or loosening fasteners.	When used as intended, the Ratcheting Driver aid in tightening or loosening fasteners.
Female Hex Driver	Remove implanted screws that have a male hex configuration.	Any surgical scenario that requires the removal of an implanted screw.	Any patient that requires the removal of an implanted screw.	The Female Hex Drivers are designed to aid in the removal of implanted screws that have a male hex configuration.	When used as intended, the Female Hex Drivers aid in the removal of implanted screws that have a male hex configuration.
Adjustable Plier	Securely grip and manipulate implants.	Any surgical scenario that requires the manipulation or gripping of a surgical rod or implant.	Any patient that requires the manipulation or gripping of a surgical rod or implant.	The Adjustable Pliers are designed to aid in manipulation or gripping of a surgical rod or implant.	When used as intended, the Adjustable Pliers aid in manipulation or gripping of a surgical rod or implant.

Product Family/Device Description	Intended Use	Indication For Use	Patient Group Target	Performance Characteristics of the Device	Expected Clinical Benefits
Screw Extractor	Remove implanted screws that have stripped or damaged heads.	Any surgical scenario that requires the removal of an implanted screw.	Any patient that requires the removal of an implanted screw.	The Screw Extractors are designed to aid in the removal of implanted screws that have stripped or damaged heads.	When used as intended, the Screw Extractors aid in the removal of implanted screws that have stripped or damaged heads.
Trephine, Broken Screw Remover, Hollow Reamer	Clear bone growth and remove a variety of damaged implanted screws.	Any surgical scenario that requires the removal of an implanted screw that may or may not be obstructed by bone-ingrowth, or may be damaged.	Any patient that requires the removal of an implanted screw that may or may not be obstructed by bone-ingrowth, or may be damaged.	Trephines, Broken Screw Removers, and Hollow Reamers are designed to aid in clearing bone growth and removing a variety of damaged implanted screws	When used as intended, the Trephines, Broken Screw Removers, and Hollow Reamers aid in clearing bone growth and removing a variety of damaged implanted screws
Custom Drivers	Remove various implant company specific screws and hardware.	Any surgical scenario that requires the removal of an implanted screw.	Any patient that requires the removal of an implanted screw.	The Custom Drivers are designed to aid in the removal of implanted screws.	When used as intended, the Custom Drivers aid in the removal of implanted screws.
Fixed Handles	Provide a torsional load to insert or remove fasteners or implants.	Any surgical scenario that requires the insertion or removal of an implanted screw.	Any patient that requires the insertion or removal of an implanted screw.	The Fixed Handles are designed to aid in the insertion or removal of fasteners or implants.	When used as intended, the Fixed Handles aid in the insertion or removal of fasteners or implants.

<b>Product Family/Device Description</b>	<b>Intended Use</b>	<b>Indication For Use</b>	<b>Patient Group Target</b>	<b>Performance Characteristics of the Device</b>	<b>Expected Clinical Benefits</b>
Counter Torque	Provide a reactionary torque to a spinal implant.	Any surgical scenario that requires a reactionary torque to be applied to a surgical implant.	Any patient with an implant that requires a reactionary torque	The Counter Torques are designed to provide a reactionary torque to an implant	When used as intended, the Counter Torques provide a reactionary torque to an implant
Extensions	Provide additional length for various ¼ square and A-O drivers used for the insertion and removal of implanted surgical devices.	Any surgical scenario that requires the insertion/removal of a surgical screw with a driver of extended length.	Any patient that requires the insertion or removal of an implanted screw.	The Extensions are designed to aid in the insertion or removal of fasteners or implants.	When used as intended, the Extensions aid in the insertion or removal of fasteners or implants.
Open Ended Wrench	Insert or remove fasteners.	Any surgical scenario that requires the insertion or removal of a fastener.	Any patient that requires the insertion or removal of a fastener.	The Open Ended Wrenches are designed to aid in the insertion or removal of fasteners.	When used as intended, the Open Ended Wrenches aid in the insertion or removal of fasteners.
Forged Rod Cutter	Shear various sized rods.	Any surgical scenario that requires a surgical rod to be cut.	Any patient that requires a surgical rod to be implanted.	The Forged Rod Cutters are designed to aid in the shearing of various sized surgical rods.	When used as intended, the Forged Rod Cutters aid in the shearing of various sized surgical rods.

## **Intended User:**

All devices listed in the table above are prescriptive and therefore to be used by qualified orthopedic surgeons trained in the respective surgical techniques. Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.

## **Contraindications**

These instruments are contraindicated in the removal of spinal hardware systems which use a pressure fit to retain its components or which twist off the threads of the locking screws.

## **Materials & Restricted Substances:**

For indication that the devices contain a restricted substance or material of animal origin see product label.

## **Possible Adverse Effects:**

None. All devices listed in the table above are used at the discretion of qualified orthopedic surgeons trained in the respective surgical techniques.

## **Precautions**

Instruments are provided non-sterile. Clean and sterilize prior to each use. Personal Protective Equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated instruments and accessories.

Neutral pH cleaners are recommended.

Prior to use, care must be taken to protect nerves, vessels and/or organs from damage that may result from the use of these instruments.

If present, safety caps and other protective packaging material must be removed from the instruments prior to the first cleaning and sterilization. Caution should be exercised while handling, cleaning, or wiping instruments and accessories with sharp cutting edges, tips, and teeth.

## **Advise:**

Devices which cannot be reprocessed by the following processes are labeled accordingly. Individual procedures accompanying these devices must be followed. In case of damage, the device should be reprocessed before returning to the manufacturer for repair.

## **Adverse Events and Complications**

### *Reporting of Serious Incidents*

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- A) the death of a patient, user, or other person,
- b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- c) a serious public health threat.

Where further information is desired, please contact your local Tecomet sales representative. For instruments produced by another legal manufacturer, reference the manufacturer's instructions for use.

## **Reprocessing Limitations**

Repeated processing has minimal effect upon manual instruments unless otherwise noted. End of life is determined by wear and damage incurred during the intended surgical use.

Non-foaming, neutral pH enzymatic and cleaning agents are recommended for processing reusable instruments and accessories.

In countries where reprocessing requirements are more stringent than those provided in this document it is the responsibility of the user/processor to comply with those prevailing laws and ordinances.

These reprocessing instructions have been validated as being capable or preparing the instruments and accessories for surgical use. It is the responsibility of the user/hospital/healthcare provider to ensure that reprocessing is performed using the appropriate equipment, materials, and that personnel have been adequately trained in order to achieve the desired result; this normally requires that the equipment and processes are validated and routinely monitored. Any deviation by the user/hospital/healthcare provider from these instructions should be evaluated for effectiveness to avoid potential adverse consequences.

### **Point of Use**

Remove excess bodily fluids and tissue with a disposable wipe. Instruments should not be allowed to dry prior to cleaning.

### **Containment / Transportation**

Follow hospital protocols when handling contaminated and/or biohazardous materials. Instrument should be cleaned within 30 minutes after use to minimize the potential of staining, damage, and drying after use.

### **Preparation for Decontamination:**

If possible, the devices must be reprocessed in a disassembled or opened state.

### **Pre-Cleaning:**

Not required.

### **Manual Cleaning Procedure**

1. Prepare a proteolytic enzyme solution according to the manufacturer's instructions.

2. Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces.
3. Soak instruments for a minimum of 10 minutes. While soaking, scrub surfaces using a nylon-bristled brush until all visible soil has been removed. Actuate moveable mechanisms. Particular attention should be given to crevices, hinged joints, box locks, instrument teeth, rough surfaces, and areas with moving components or springs. Lumens, blind holes, and cannulas should be cleaned using a snug fitting round nylon bristle brush. Insert the snug fitting round brush into the lumen, blind hole, or cannula with a twisting motion while pushing in and out multiple times.

***NOTE: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.***

4. Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas, and other difficult to access areas.
5. Prepare an ultrasonic cleaning bath with detergent and de-gas according to the manufacturer's recommendations. Completely submerge instruments in the cleaning solution and gently shake them to remove any trapped bubbles. Lumens, blind holes, and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces. Sonically clean the instruments at the time, temperature, and frequency recommended by the equipment manufacturer and

optimal for the detergent used. A minimum of ten (10) minutes is recommended.

**NOTES:**

- *Separate stainless steel instruments from other metal instruments during ultrasonic cleaning to avoid electrolysis.*
  - *Fully open hinged instruments.*
  - *Use wire mesh baskets or tray designed for ultrasonic cleaners.*
  - *Regular monitoring of sonic cleaning performance by means of an ultrasonic activity detector, aluminum foil test, TOSI™ or SonoCheck™ is recommended.*
6. Remove the instruments from the ultrasonic bath and rinse in purified water for a minimum of one (1) minute or until there is no sign of residue detergent or biologic soil. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas, and other difficult to access areas.
  7. Dry instruments with a clean, absorbent, non-shedding lint free cloth. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas, and difficult to access areas.

**Combination Manual/Automated Cleaning Procedure**

1. Prepare a proteolytic enzyme solution according to the manufacturer's instructions.
2. Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes, and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces.
3. Soak instruments for a minimum of ten (10) minutes. Scrub surfaces using a soft nylon-bristled brush until all visible soil has been removed. Actuate moveable mechanisms. Particular attention should be given to crevices, hinged joints, box locks,

instrument teeth, roughened surfaces, and areas with moving components or springs. Lumens, blind holes and cannulas should be cleaned using a snug fitting round nylon bristle brush. Insert the snug fitting round nylon bristle brush into the lumen, blind hole or cannula with a twisting motion while pushing in and out multiple times.

**Note: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.**

4. Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas, and other difficult to access areas.
5. Place instruments in a suitable washer/disinfector. Follow the washer/disinfector manufacturer's instructions for loading the instruments for maximum cleaning exposure; e.g. open all instruments, place concave instrument on their side or upside down, use baskets and trays designed for washers, place heavier instruments on the bottom of trays and baskets. If the washer/disinfector is equipped with special racks (e.g. for cannulated instruments) use them according to the manufacturer's instructions.

6. Process instruments using a standard washer/disinfector instrument cycle according to the manufacturer's instructions. The following minimum wash cycle parameters are recommended:

Cycle	Description
1	Pre-wash • Cold Softened Tap Water • 2 minutes
2	Enzyme Spray & Soak • Hot Softened Tap Water • 1 minute
3	Rinse • Cold Softened Tap Water • 1 minute
4	Detergent Wash • Hot Tap Water (64-66°C/146-150°F) • 2 minutes
5	Rinse • Hot Purified Water (64-66°C/146-150°F) • 1 minute
6	Hot Air Dry (116°C/240°F) • 7 – 30 minutes

**Notes:**

- The washer/disinfector manufacturer's instructions should be followed.
- A washer/disinfector with demonstrated efficacy (e.g. FDA approved, validated to ISO 15883) should be used.
- Dry time is shown as a range because it is dependent upon the load size placed into the washer/disinfector.
- Many manufacturers pre-program their washer/disinfectors with standard cycles and they may include a thermal low-level disinfection rinse after the detergent wash. The thermal disinfection cycle should be performed to achieve a minimum value A0=600(e.g. 90°C/194°F for 1 minute according to ISO 15883-1) and is compatible with instruments.
- If a lubrication cycle is available that applies to water-soluble lubricant such as Preserve®, Instrument Milk or equivalent material intended for medical device application it is acceptable to use on Instruments unless otherwise indicated.

## **Disinfection**

Disinfection is only acceptable as an adjunct to full sterilization for reusable surgical instruments. See sterilization section below.

Low level disinfection may be used as part of a washer/disinfector cycle but the devices must also be sterilized before use.

## **Drying:**

Dry instruments with a clean, absorbent, non-shedding lint free cloth. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas, and difficult to access areas.

## **Inspection / Functional Testing**

1. Carefully inspect each device to ensure that all visible blood and soil has been removed.
2. Visually inspect for damage and/or wear.
3. Check the action of moving parts (such as hinges and box-locks) to ensure smooth operation throughout the intended range of motion.
4. Check instruments with long slender features (particularly rotating instruments) for distortion.
5. Where instruments form part of a larger assembly, check that the devices assemble readily with mating components.
6. Inspect edges of the case, tray and cassette to ensure there aren't any sharp edges which might cut sterilization wrap. Inspect lid to assure there isn't any cracking, bending or warping. Verify lid continues to fit case tightly.

***\*Note: If damage or wear is noted that may compromise the function of the instrument, do not use the instrument, and notify the appropriate responsible person.***

## **Sterilization**

Instruments and the instrument kit must be cleaned prior to sterilization. Instruments must be arranged per the specified tray configuration. Additional instruments not included in the specified tray configuration shall not be added to the tray for sterilization.

With instruments secured within, the kit must be wrapped using CSR (sterilization wrap) following the AAMI double wrap method (ANSI/AAMI ST46), or equivalent. The trays and lids of the kit must be properly assembled prior to sterilization. It is not recommended to stack assembled kits during sterilization. Do not place kit on side or upside down. The use of nonabsorbent tray liners can cause condensation to pool.

Single devices may be packaged in an approved (e.g. FDA cleared or ISO 11607 compliant) medical grade sterilization pouch or wrap. Care should be used when packaging so that the pouch or wrap is not torn. Devices should be wrapped using the double wrap or equivalent method (ref: AAMI ST79, AORN Guidelines).

Reusable wraps are not recommended.

Moist heat/steam sterilization is the recommended method for the instruments. Use of an approved chemical indicator (class 5) or chemical emulator (class 6) within each sterilization load is recommended. Always consult and follow the sterilizer manufacturer instructions for load configuration and equipment operation. Sterilizing equipment should have demonstrated efficacy (e.g. FDA clearance, EN 13060 or EN 285 compliance). Additionally the manufacturer's recommendations for installation, validation, and maintenance should be followed.

Validated exposure time and temperature to achieve a  $10^{-6}$  sterility assurance level (SAL) is listed below. Local or national specification should be followed where steam sterilization requirements are stricter or more conservative than those listed below.

**Method – Pre-Vacuum Steam Sterilization – for 2nd Generation Evolution Systems (Part Numbers ES3XXXX and ES4XXXX ONLY. For example ES30001 or ES40001)**

Set parameters at 4 minutes exposure at a temperature of 132°C (270°F) with four (4) vacuum pulses prior to actual exposure. A dry time of 30 minutes is recommended.

**Method – Gravity Displacement Steam Sterilization – for 2nd Generation Evolution Systems (Part Numbers ES3XXXX or ES4XXXX ONLY. For example ES30001 or ES40001)**

Expose instruments for 15 minutes at a temperature of 132° C (270° F). A dry time of 30 minutes is recommended.

**Method – Pre-Vacuum Steam Sterilization-for 1st Generation Evolution Systems (Any other Part Number not listed above.)**

Set parameters at 10 minutes exposure at a temperature of 132°C (270°F) with four (4) vacuum pulses prior to actual exposure.

**Method – Gravity Displacement Steam Sterilization- for 1st Generation Evolution Systems (Any other Part Number not listed above)**

Expose instruments for 30 minutes at a temperature of 132° C (270° F).

Handle sterilized kit following normal operating room procedures. Allow case to cool thoroughly for a minimum of 15 minutes after sterilization cycle is complete.

**Maintenance**

After cleaning and before sterilization, instruments may be lubricated with a water-soluble lubricant such as Preserve®, Instrument Milk or equivalent material intended for surgical application. Lubricate hinges, threads and other moving part. Always follow the lubricant manufacturer's instructions for dilution, shelf life and application method.

## **Storage**

Store sterile, packaged instruments in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes. **Note:** Inspect every package before use to ensure that the sterile barrier (e.g. wrap, pouch or filter) is not torn, perforated, shows signs of moisture, or appears to be tampered with. If any of those conditions are present then the contents are considered non-sterile and should be reprocessed through cleaning, packaging, and sterilization.

## **Instructions for use**

**Cautions:** The user of this instrument set acknowledges that implant manufacturers are constantly improving their designs, which could lead to a fit concern with custom drivers.

Do not use the driver(s) if the instrument does not properly fit into the implanted spinal hardware.

Prior to inserting the driver into the implant, verify that all foreign material is removed from the driving surface. Failure to do so could result in the damaging of implants and/or removal instruments.

The user of the instrument set should ensure that the instrumentation has no more than 60 minutes of direct exposure to the patient. The user of this set should ensure that no instrumentation/parts of instrumentation be left in the patient during and following the surgical procedure.

Inspect instruments for damage or wear prior to use. Failure to do so could result in the damaging of implants or removal instruments.

**\*Note: If damage or wear is noted that may compromise the function of the instrument, do not use the instrument and notify the appropriate responsible person.**

The following tables identify the torque ratings of the male hex and male hexalobe drivers. The drivers should not be subjected to torque loads higher than the identified ratings in order to ensure safe use.

### **Male Hex Driver Torque Ratings**


<b>Size</b>	<b>Torque Rating (in.-lbs)</b>
1.5 mm	7
2.0 mm	10
2.5 mm	29
7/64"	29
3.0 mm	47
3.2 mm (1/8")	47
3.5 mm	71
5/32"	100
4.5 mm	138
3/16"	138
5.0 mm	224
7/32"	224
1/4"	224
7.0 mm	224

### **Male Hexalobe Driver Torque Ratings**

<b>Size</b>	<b>Torque Rating (in.-lbs)</b>
T08	10
T10	17
T15	29
T20	47
T25	71
T27	100
T30	138
T40	240

## BROKEN SCREW REMOVAL


### (Combined trephine and Bone Screw Remover)

1. Measure the outside diameter of broken screw and select trephine/Broken Screw Remover that matches screw size.
2. Place trephine in handle. Turn counter-clockwise around broken screw. The trephine should clear around the screw as well as bite into screw.
3. Continue turning in counter-clockwise direction until screw is completely removed.

***Note: There may be bone in-growth on the screw preventing removal of the screw. Depending on the surgical scenario and condition of implanted screws, the broken screw remover may become damaged after normal usage. Inspect the broken screw remover following usage to determine if any damage has occurred. Discontinue use if damage is present.***

## BROKEN SCREW REMOVAL (Separate trephine and Bone Screw Remover)

1. Measure the outside diameter of the screw and select trephine that matches screw size.
2. Place trephine in handle of choice provided in the removal kit. Position trephine over broken screw. Note: The Trepbine should clear around the screw. Approximately 5mm of screw should be exposed to properly seat the Trepbine. To expose the broken screw, use a clockwise and counter-clockwise motion while applying downward pressure. Advance the trephine in this manner until the trephine reaches the broken screw. Note: During this process the teeth of the trephine may become clogged with debris. Remove trephine and clean debris from tool.
3. Select the bone screw remover that matches the diameter of the bone screw to be removed. Place bone screw remover into handle and position bone screw remover over the exposed screw. Use a counter-clockwise motion while applying downward pressure. The bone screw remover will bite into the screw. Continue turning in a counter-clockwise direction until the screw is completely removed.

***Note: There may be bone in-growth on the screw preventing removal of the screw. Depending on the surgical scenario and condition of implanted screws, the broken screw remover may become damaged after normal usage. Inspect the broken screw remover following usage to determine if any damage has occurred. Discontinue use if damage is present.***

The male hex drivers, female hex drivers, male hexalobe drivers, and custom drivers are to be used with the appropriate  $\frac{1}{4}$ " square or AO ratcheting or fixed driver to remove implanted screws. The driver type should be selected based on the appropriate connection style of the implanted screw. Insert the appropriate driver tip into the mating features of the implanted screw and rotate the driver in the counter-clockwise direction to remove the screw.

The screw extractors can be used to remove stripped implanted screws. The screw extractors are to be used with the appropriate  $\frac{1}{4}$ " square or AO ratcheting or fixed driver. To remove a stripped screw, rotate the screw extractor in the counter-clockwise direction on the stripped screw head surface.

The ratcheting driver is to be used with the appropriate bit driver needed for screw removal. The bit driver can be inserted into the ratcheting driver by compressing the adaptor collar and inserting the appropriate  $\frac{1}{4}$ " square or AO driver. Ensure that the driver is fully engaged prior to use. The ratchet can operate in a forward or reverse ratcheting direction, or in a locked setting (fixed driver). The direction of the ratcheting driver can be changed by rotating the selector ring in the clockwise or counter-clockwise direction. The bit driver can be removed from the ratcheting driver by compressing the adaptor collar.

The fixed driver is to be used with the appropriate bit driver needed for screw removal. The bit driver can be inserted into the fixed driver by compressing the adaptor collar and inserting the appropriate  $\frac{1}{4}$ " square or AO driver. The bit driver can be removed from the fixed driver by compressing the adaptor collar.

The extensions are to be used with a fixed or ratcheting handle to provide additional length to the bit drivers used to remove implanted screws. The extensions can be inserted into the adaptors of the ratcheting or fixed drivers. Bit drivers can be inserted into the extensions by either compressing the adaptor collar and inserting (AO extensions) or by inserting the male bit driver connection end into the female connection end of the extension (1/4" sq extensions). The bit drivers can be removed from the extensions by either compressing the adaptor collar (AO extensions) or pulling the bit driver out of engagement (1/4" sq extensions).

The counter torque wrenches are used to apply a counter torque or opposite torque to a spinal rod during surgery. The counter torque can be inserted into the appropriate fixed handle. The tip of the counter torque with the full radius slot can be placed over the rod to provide stability as needed.

The adjustable pliers can be used to grip spinal rods during the screw extraction process. The pliers can lock based on the position of the screw knob on one of the handles.

The open ended wrenches are used to remove fasteners. The open ended wrenches can remove various sized fasteners based on the open ended wrench size. Engage the wrench onto the mating features of the implanted screw and rotate the wrench in the counter-clockwise direction to remove the screw.

The tri-lobe adapter can be used with a handle with a female tri-lobe adapter or a 3-jaw chuck adapter. An appropriate sized bit driver can be inserted into the female square end of the tri-lobe adapter in order to remove an implanted screw.

The rod cutter can be used to shear rods in various sizes, as dictated by the surgical procedure and surgeon. To cut a rod, open the rod cutter handles by pulling away from each other, insert a rod to be lined up with the rod cutter teeth, and apply force to close the rod cutter handles by pushing the handles towards each other.

**Disposal:**

Devices are reusable; but may be disposed of at the end of the device life by means of standard hospital disposal procedures.

**Additional Advices:**

It is the duty of the user to ensure that the reprocessing processes are followed; that resources and materials are available to capable personnel; and that hospital protocols are followed. State of the art and often national law requiring these processes and included resources to be validated and maintained properly.

## Symbols Used on Labeling<sup>1</sup>:

	Caution
	Non-Sterile
	Sterile
	U.S Federal law restricts this device to sale by or on the order of a physician
	CE Mark <sup>1</sup>
	CE Mark with Notified Body # <sup>1</sup>
	Authorized Representative in the European Community
	Manufacturer
	Date of Manufacture
	Lot Number
	Catalog Number
	Consult Instructions for Use
	Medical Device
	Packaging Unit
	Country of Manufacture
	Distributor
	Swiss Authorized Representative

Importer

Unique Device Identifier

**Additional Labeling:**

**"MANUAL USE  
ONLY"**

Device should not be attached to a power source and is only intended to be handled manually.

**"REMOVE FOR  
CLEANING &  
STERILIZATION"**

Device should be disassembled prior to cleaning and sterilization.

<sup>1</sup>Refer to the labelling for CE information

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Symmetry Medical Manufacturing, Inc.
S.A.S