TECOMET

Manufactured By:
TECOMET
 Tecomet, Inc.
 5307 95th Avenue
 Kenosha, WI 53144 USA
 Toll Free 855-925-4500
 262-925-4500

www.tecomet.com

EC REP European Representative
 Meditec Source GmbH & Co. KG
 Sattlerstraße 19 78532
 Tuttlingen, Germany

ENGLISH	Reprocessing Surgical Instruments	やまと	手術器具の再処理
БЪЛГАРСКИ	Повторна обработка на хирургически инструменти	LATVIĒSU	Kirurgisko instrumentu atkārtota apstrāde
HRVATSKI	Ponovna obrada kirurških instrumenata	LIETUVIŲ	Chirurginių instrumentų pakartotinis apdorojimas
CZECH	Obnova chirurgických nástrojů	NORSK	Repressering av kirurgiske instrumenter
DANSK	Oparbejdning af kirurgiske instrumenter	POLSKI	Regeneracja narzędzi chirurgicznych
NEDERLANDS	Herverwerking van chirurgische instrumenten	PORTUGUÊS	Reprocessamento de instrumentos cirúrgicos
EESTI	Kirurgiliste instrumentide taastõtlemine	ROMÂNĂ	Reprelucrarea instrumentelor chirurgicale
SUOMI	Kirurgisten instrumenttien uudelleenkäsitely	SLOVENSKY	Renovanie chirurgických nástrojov Slovensky
FRANÇAIS	Retraitement des instruments chirurgicaux	SLOVENŠČINA	Kirurški instrumenti za ponovno obdelavo
DEUTSCH	Aufbereitung chirurgischer Instrumente	ESPAÑOL	Reprocesamiento de instrumentos quirúrgicos
ΕΛΛΗΝΙΚΗ	Επανεξεργασία χειρουργικών εργαλείων	SVENSKA	Reprocessing kirurgiska instrument
MAGYAR	Sebészeti műszerek újrafeldolgozása	TÜRKÇE	Cerrahi Aletleri Tekrar İşleme
ITALIANO	Ricondizionamento degli strumenti chirurgici		

Reprocessing Surgical Instruments

Issue Date: 2022/02/24

Product

Surgical instruments

Manufacturer

Tecomet, Inc.
5307 95th Avenue
Kenosha, WI 53144 USA

European Representative:

Meditec Source GmbH & Co. KG
Sattlerstraße 19 78532
Tuttlingen, Germany

These instructions are in accordance with ISO17664 and AAMI ST81. They apply to reusable surgical instruments and accessories (provided as non-sterile) supplied by Tecomet and intended for reprocessing in a health care facility setting. All instruments and accessories may be safely and effectively reprocessed using the manual or combination manual/automated cleaning instructions and sterilization parameters provided in this document **UNLESS otherwise noted in instructions accompanying a specific instrument.**

Intended Use:

Table 1: Product family information for non-measuring devices

Product Family/ Device Description	Intended Use	Indication For Use	Patient Group Target	Performance Characteristics of the Device	Expected Clinical Benefits
Male Hexalobe Driver	Remove implanted screws that have a female Hexalobe configuration.	Any surgical scenario that requires the removal of an implanted screw.	Any patient that requires the removal of an implanted screw.	The Male Hexalobe Drivers are designed to aid in the removal of implanted screws that have a female Hexalobe configuration.	When used as intended, the Male Hexalobe Drivers, aid in the removal of implanted screws that have a female Hexalobe configuration.
Male Hex Driver	Remove implanted screws that have a female hex configuration.	Any surgical scenario that requires the removal of an implanted screw.	Any patient that requires the removal of an implanted screw.	The Male Hex Drivers are designed to aid in the removal of implanted screws that have a female hex configuration.	When used as intended, the Male Hex Drivers, aid in the removal of implanted screws that have a female hex configuration.
Ratcheting Driver	Increase efficiency when tightening or loosening fasteners.	Any surgical scenario where a clinician desires to increase efficiency when tightening or loosening fasteners.	Any patient that requires tightening or loosening of fasteners.	The Ratcheting Driver are designed to aid in tightening or loosening fasteners.	When used as intended, the Ratcheting Driver aid in tightening or loosening fasteners.
Female Hex Driver	Remove implanted screws that have a male hex configuration.	Any surgical scenario that requires the removal of an implanted screw.	Any patient that requires the removal of an implanted screw.	The Female Hex Drivers are designed to aid in the removal of implanted screws that have a male hex configuration.	When used as intended, the Female Hex Drivers aid in the removal of implanted screws that have a male hex configuration.
Adjustable Plier	Securely grip and manipulate implants.	Any surgical scenario that requires the manipulation or gripping of a surgical rod or implant.	Any patient that requires the manipulation or gripping of a surgical rod or implant.	The Adjustable Pliers are designed to aid in manipulation or gripping of a surgical rod or implant.	When used as intended, the Adjustable Pliers aid in manipulation or gripping of a surgical rod or implant.
Screw Extractor	Remove implanted screws that have stripped or damaged heads.	Any surgical scenario that requires the removal of an implanted screw.	Any patient that requires the removal of an implanted screw.	The Screw Extractors are designed to aid in the removal of implanted screws that have stripped or damaged heads.	When used as intended, the Screw Extractors aid in the removal of implanted screws that have stripped or damaged heads.
Trephine, Broken Screw Remover, Hollow Reamer	Clear bone growth and remove a variety of damaged implanted screws.	Any surgical scenario that requires the removal of an implanted screw that may or may not be obstructed by bone-ingrowth, or may be damaged.	Any patient that requires the removal of an implanted screw that may or may not be obstructed by bone-ingrowth, or may be damaged.	Trephines, Broken Screw Removers, and Hollow Reamers are designed to aid in clearing bone growth and removing a variety of damaged implanted screws	When used as intended, the Trephines, Broken Screw Removers, and Hollow Reamers aid in clearing bone growth and removing a variety of damaged implanted screws

Product Family/ Device Description	Intended Use	Indication For Use	Patient Group Target	Performance Characteristics of the Device	Expected Clinical Benefits
Custom Drivers	Remove various implant company specific screws and hardware.	Any surgical scenario that requires the removal of an implanted screw.	Any patient that requires the removal of an implanted screw.	The Custom Drivers are designed to aid in the removal of implanted screws.	When used as intended, the Custom Drivers aid in the removal of implanted screws.
Fixed Handles	Provide a torsional load to insert or remove fasteners or implants.	Any surgical scenario that requires the insertion or removal of an implanted screw.	Any patient that requires the insertion or removal of an implanted screw.	The Fixed Handles are designed to aid in the insertion or removal of fasteners or implants.	When used as intended, the Fixed Handles aid in the insertion or removal of fasteners or implants.
Counter Torque	Provide a reactionary torque to a spinal implant.	Any surgical scenario that requires a reactionary torque to be applied to a surgical implant.	Any patient with an implant that requires a reactionary torque	The Counter Torques are designed to provide a reactionary torque to an implant	When used as intended, the Counter Torques provide a reactionary torque to an implant
Extensions	Provide additional length for various 1/4 square and A-O drivers used for the insertion and removal of implanted surgical devices.	Any surgical scenario that requires the insertion/removal of a surgical screw with a driver of extended length.	Any patient that requires the insertion or removal of an implanted screw.	The Extensions are designed to aid in the insertion or removal of fasteners or implants.	When used as intended, the Extensions aid in the insertion or removal of fasteners or implants.
Open Ended Wrench	Insert or remove fasteners.	Any surgical scenario that requires the insertion or removal of a fastener.	Any patient that requires the insertion or removal of a fastener.	The Open Ended Wrenches are designed to aid in the insertion or removal of fasteners.	When used as intended, the Open Ended Wrenches aid in the insertion or removal of fasteners.
Forged Rod Cutter	Shear various sized rods.	Any surgical scenario that requires a surgical rod to be cut.	Any patient that requires a surgical rod to be implanted.	The Forged Rod Cutters are designed to aid in the shearing of various sized surgical rods.	When used as intended, the Forged Rod Cutters aid in the shearing of various sized surgical rods.
Bi-Directional Palm Torque Limiting Instruments	Limit the applied amount of torque in the forward and/or reverse direction to a pre-determined calibrated setting.	Any surgical scenario where a clinician desires to apply a specific torsional load to a fastener or implant in either the clockwise or counter-clockwise direction.	Any patient that requires a fastener or implant set to a specific torsional load.	The Bi-Directional Palm Torque Limiting Instruments are designed to apply a specific torsional load to a fastener or implant	When used as intended, the Bi-Directional Palm Torque Limiting Instruments apply a specific torsional load to a fastener or implant
Silicone Overmolded Handles with Lasermarking Options	Provide an area for the surgeon to grip and apply a load to an instrument, accessory and/or implant as a component of a finished instrument.	Any surgical scenario where a clinician desires to apply a load to an instrument, accessory and/or implant.	Any patient that requires an instrument, accessory and/or implant which needs load applied	The Silicone Overmolded Handles with Lasermarking Options are designed to provide an area for the surgeon to grip and apply a load as a component of a finished instrument.	When used as intended, the Silicone Overmolded Handles with Lasermarking Options provide an area for the surgeon to grip and apply a load as a component of a finished instrument.

Product Family/ Device Description	Intended Use	Indication For Use	Patient Group Target	Performance Characteristics of the Device	Expected Clinical Benefits
Sounders	Aid in the detection and differentiation of varying tissue types.	Any surgical scenario where a clinician desires to differentiate tissues.	Any patient undergoing a surgery where a clinician desires to differentiate tissues.	The Sounders are designed to aid in the detection and differentiation of varying tissue types.	When used as intended, the Sounders aid in the detection and differentiation of varying tissue types.
Awls	Pierce holes into bone to initiate pilot holes for surgical screws or fasteners.	Any surgical scenario where a clinician must pierce holes into bone.	Any patient undergoing a surgery where a clinician must pierce holes into bone.	The Awls are designed to aid in piercing holes into bone to initiate pilot holes for surgical screws or fasteners.	When used as intended, the Awls aid in piercing holes into bone to initiate pilot holes for surgical screws or fasteners.
Curettes	Allow for incremental removal of various tissues.	Any surgical scenario where a clinician desires to scrape at and remove tissues.	Any patient undergoing a surgery where a clinician desires to scrape at and remove tissues.	The Curettes are designed to aid in incremental removal of various tissues.	When used as intended, the Curettes aid in incremental removal of various tissues.
Rod Cutters	Shear various sized rods based on the clinician and surgeon's needs.	Any surgical scenario that requires a surgical rod to be cut.	Any patient that requires a surgical rod to be implanted.	The Rod Cutters are designed to aid in the shearing of various sized surgical rods.	When used as intended, the Rod Cutters aid in the shearing of various sized surgical rods.
Rod Benders	Bend various sized rods based on the clinician and surgeon's needs.	Any surgical scenario that requires a surgical rod to be bent.	Any patient that requires a surgical rod to be implanted.	The Rod Benders are designed to aid in the bending of various sized surgical rods.	When used as intended, the Rod Benders aid in the bending of various sized surgical rods.
Compressors and Distractors	Apply compression or distraction at various vertebrae levels with various sized spinal rods prior to final tightening of the set screw.	Any surgical scenario where a clinician desires to increase efficiency when compressing, distracting, or spreading various materials.	Any patient undergoing a surgery where a clinician desires to increase efficiency when compressing, distracting, or spreading various materials.	The Compressors and Distractors are designed to aid in applying compression or distraction at various vertebrae levels with various sized spinal rods prior to final tightening of the set screw.	When used as intended, the Compressors and Distractors aid in applying compression or distraction at various vertebrae levels with various sized spinal rods prior to final tightening of the set screw.
Plate Benders	Bend various sized plates.	Any surgical scenario that requires a surgical plate to be bent.	Any patient that requires a surgical plate to be implanted.	The Plate Benders are designed to aid in the bending of various sized surgical plates.	When used as intended, the Plate Benders aid in the bending of various sized surgical plates.
Wire Cutters	Shear various sized wires.	Any surgical scenario that requires a surgical wire to be cut.	Any patient that requires a surgical wire to be used or implanted.	The Wire Cutters are designed to aid in the shearing of various sized surgical wire.	When used as intended, the Wire Cutters aid in the shearing of various sized surgical wire.

Product Family/ Device Description	Intended Use	Indication For Use	Patient Group Target	Performance Characteristics of the Device	Expected Clinical Benefits
Trials and Templates	Simulate an implant that will be used during surgery. The size and shape will be based on the clinician and surgeon's needs.	Any surgical scenario that requires a representation of an implant as based on the clinician and surgeon's needs.	Any patient that requires an implant.	The Trials and Templates are designed to simulate an implant that will be used during surgery.	When used as intended, the Trials and Templates simulate an implant that will be used during surgery.
Mallets	Create impact force.	Any surgical scenario where a clinician desires to create an impact force.	Any patient undergoing a surgery where a clinician desires to create an impact force.	The Mallets are designed to aid in creating an impact force.	When used as intended, the Mallets aid in creating an impact force
Rasps	File or coarsely shape bone or hard tissue.	Any surgical scenario where a clinician intends to remove, file, or shape bone or hard tissue.	Any patient undergoing a surgery where a clinician intends to remove, file, or shape bone or hard tissue.	The Rasps are designed to aid in filing or coarsely shaping bone or hard tissue.	When used as intended, the Rasps aid in filing or coarsely shaping bone or hard tissue.
Bone Tamps	Manipulate, support, and impact bone.	Any surgical scenario where a clinician desires to transfer axial load.	Any patient undergoing a surgery where a clinician desires to transfer axial load.	The Bone Tamps are designed to aid in manipulating, supporting, and impacting bone.	When used as intended, the Bone Tamps aid in manipulating, supporting, and impacting bone.
Forceps	Compress or grasp bone, tissue, or hardware.	Any surgical scenario where a clinician desires to compress or grasp bone, tissue, hardware, or instrumentation.	Any patient undergoing a surgery where a clinician desires to compress or grasp bone, tissue, hardware, or instrumentation.	The Forceps are designed to aid in compressing or grasping bone, tissue, or hardware.	When used as intended, the Forceps aid in compressing or grasping bone, tissue, or hardware.
Rongeurs	Remove bone and/or tissue.	Any surgical scenario that requires the removal of bone and/or tissue.	Any patient undergoing a surgery where a clinician requires the removal of bone and/or tissue.	The Rongeurs are designed to aid in removing bone and/or tissue	When used as intended, the Rongeurs aid in removing bone and/or tissue
Manipulative Instruments	Maneuver and hold tissue and/or bone.	Any surgical scenario where a clinician desires to displace and hold tissues and bones that obstruct the surgical site.	Any patient undergoing a surgery where a clinician desires to displace and hold tissues and bones that obstruct the surgical site.	The Manipulative Instruments are designed to aid in maneuvering and holding tissue and/or bone.	When used as intended, the Manipulative aid in maneuvering and holding tissue and/or bone.
Scissors and Shears	Cut through tissue and bandages.	Any surgical scenario where a clinician desires to cut bandages and tissue.	Any patient undergoing a surgery where a clinician desires to cut bandages and tissue.	The Scissors and Shears are designed to aid in cutting through tissue and bandages.	When used as intended, the Scissors and Shears aid in cutting through tissue and bandages.
Osteotomes and Gouges	Cut, dissect, remove, or prepare bone.	Any surgical scenario where a clinician desires to cut, dissect, or prepare bone.	Any patient undergoing a surgery where a clinician desires to cut, dissect, or prepare bone.	The Osteotomes and Gouges are designed to aid in cutting, dissecting, removing, or preparing bone.	When used as intended, the Osteotomes and Gouges aid in cutting, dissecting, removing, or preparing bone.

Product Family/ Device Description	Intended Use	Indication For Use	Patient Group Target	Performance Characteristics of the Device	Expected Clinical Benefits
Drill Guides	Aid and control drilling and/or insertion of bone screws.	Any surgical scenario where a clinician desires to create holes or cuts in bone.	Any patient undergoing a surgery where a clinician desires to create holes or cuts in bone.	The Drill Guides are designed to aid and control drilling and/ or insertion of bone screws.	When used as intended, the Drill Guides aid and control drilling and/ or insertion of bone screws.
Countersinks, Drill Bits, and Cutters	Remove bone through a rotational cut.	Any surgical scenario where a clinician desires to create holes or cuts in bone.	Any patient undergoing a surgery where a clinician desires to create holes or cuts in bone.	The Countersinks, Drill Bits, and Cutters are designed to aid in removing bone through a rotational cut.	When used as intended, the Countersinks, Drill Bits, and Cutters aid in removing bone through a rotational cut.
Holders and Sleeves	Aid in the insertion and removal of screws.	Any surgical scenario where a clinician desires to insert or remove a screw in a surgical setting.	Any patient undergoing a surgery where a clinician desires to insert or remove a screw in a surgical setting.	The Holders and Sleeves are designed to aid in the insertion and removal of screws.	When used as intended, the Holders and Sleeves aid in the insertion and removal of screws.
Legacy Torque Drivers	Apply a specific torsional load to a fastener or implant in the clockwise direction.	Any surgical scenario where a clinician desires to apply a specific torsional load to a fastener or implant in a clockwise direction.	Any patient that requires a fastener or implant set to a specific torsional load.	The Legacy Torque Drivers are designed to apply a specific torsional load to a fastener or implant	When used as intended, the Legacy Torque Drivers apply a specific torsional load to a fastener or implant
Rod Pushers	Push or manipulate rods.	Any surgical scenario where a clinician desires to push or manipulate rods.	Any patient undergoing a surgery where a clinician desires to push or manipulate rods.	The Rod Pushers are designed to aid in pushing or manipulating rods	When used as intended, the Rod Pushers aid in pushing or manipulating rods
Guidewire Holders	Securely grip and hold wires and/or guide wires.	Any scenario where a surgeon desires to securely grip and hold wires and/or guide wires into a surgical site.	Any patient that requires a surgical wire to be used or implanted.	The Guidewire Holders are designed to aid in securely gripping and holding wires and/or guide wires.	When used as intended, the Guidewire Holders aid in securely gripping and holding wires and/or guide wires.
Angled Drivers	Translate or provide a torsional load to insert or remove fasteners or implants using an angled head feature	Any surgical scenario where a clinician desires to translate or provide a torsional load to insert or remove fasteners or implants using an angled head feature.	Any patient that requires the insertion or removal of an implanted screw.	The Angled Drivers are designed to aid in translating or providing a torsional load to insert or remove fasteners or implants	When used as intended, the Angled Drivers aid in translating or providing a torsional load to insert or remove fasteners or implants
Deflection Beam Torque Instruments	Allow a surgeon to apply a specific torsional load to a fastener or implant.	Any surgical scenario where a clinician desires to apply a specific torsional load to a fastener.	Any patient that requires a fastener or implant set to a specific torsional load.	The Deflection Beam Torque Instruments are designed to apply a specific torsional load to a fastener or implant	When used as intended, the Deflection Beam Torque Instruments apply a specific torsional load to a fastener or implant

Table 2: Product family information for measuring devices including claimed measuring degree of accuracy

Product Family/Device Description	Intended Use	Indication For Use	Measuring Degree of Accuracy	Patient Group Target	Performance Characteristics of the Device	Expected Clinical Benefits
Bone Probes	Measure the depth of anatomical features and/or act as a blunt dissector.	Any surgical scenario where a clinician desires to measure the depth of anatomical features or dissect tissues.	± 0.015 in	Any patient undergoing a surgery where a clinician desires to measure the depth of anatomical features and/or dissect tissues.	The Bone Probes are designed to aid in measuring the depth of anatomical features and/or acting as a blunt dissector	When used as intended, the Bone Probes aid in measuring the depth of anatomical features and/or acting as a blunt dissector
Depth Gages	Measure depth of holes or grooves.	Any surgical scenario where a clinician desires to measure the depth of holes or grooves.	± 0.5 mm	Any patient undergoing a surgery where a clinician desires to measure the depth of anatomical features and/or dissect tissues.	The Bone Probes are designed to aid in measuring the depth of anatomical features and/or acting as a blunt dissector	When used as intended, the Bone Probes aid in measuring the depth of anatomical features and/or acting as a blunt dissector

Contraindications:

There are no known contraindication at this time for the devices listed in tables 1 and 2.

Materials & Restricted Substances:

For indication that the device contains a restricted substance or material of animal origin see product label.

Intended User:

All devices listed in tables 1 and 2 are prescriptive and therefore to be used by qualified orthopedic surgeons trained in the respective surgical techniques. **Caution:** Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Possible Adverse Effects:

None. All devices listed in tables 1 and 2 are used at the discretion of qualified orthopedic surgeons trained in the respective surgical techniques.

Precautions:

Instrument is provided non-sterile. Clean and sterilize prior to each use. Personal Protective Equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated instruments and accessories.

Neutral pH cleaners are recommended.

Prior to use, care must be taken to protect nerves, vessels, and/or organs from damage that may result from the use of these instruments.

If present, safety caps and other protective packaging material must be removed from the instruments prior to the first cleaning and sterilization.

Caution should be exercised while handling, cleaning, or wiping instruments and accessories with sharp cutting edges, tips, and teeth.

Advise:

Devices which cannot be reprocessed by the following processes are labeled accordingly. Individual procedures accompanying these devices must be followed. In case of damage, the device should be reprocessed before returning to the manufacturer for repair.

Adverse Events and Complications

Reporting of Serious Incidents

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Serious incident means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- a) the death of a patient, user, or other person,
- b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- c) a serious public health threat.

Where further information is desired, please contact your local Tecomet sales representative. For instruments produced by another legal manufacturer, reference the manufacturer's instructions for use.

Reprocessing Limitations:

Repeated processing according to these instructions has minimal effect upon metal reusable instruments and accessories unless otherwise noted. End of life for stainless steel or other metal surgical instruments is generally determined by wear and damage incurred during the intended surgical use.

Non-foaming, neutral pH enzymatic and cleaning agents are recommended for processing reusable instruments and accessories.

In countries where reprocessing requirements are more stringent than those provided in this document it is the responsibility of the user/processor to comply with those prevailing laws and ordinances.

These reprocessing instructions have been validated as being capable of preparing the instruments and accessories for surgical use. It is the responsibility of the user/hospital/health care provider to ensure that reprocessing is performed using the appropriate equipment, materials and that personnel have been adequately trained in order to achieve the desired result; this normally requires that the equipment and processes are validated and routinely monitored. Any deviation by the user/hospital/health care provider from these instructions should be evaluated for effectiveness to avoid potential adverse consequences.

Regarding Legacy Torque Limiting Drivers the re-calibration interval of this product is based on its application and use in the field. It is the OEM/user/hospital/health care provider's determination as to what an appropriate re-calibration cycle is. In general, Tecomet estimates a re-calibration period as 6 month period or every 100 autoclave cycles, whichever occurs first. This is only an estimate and the OEM/user/hospital/health care provider must determine the appropriate recalibration cycle based on the specific circumstances, applications, and scenarios where a torque limiting driver would be utilized.

Preparation at the Point of Use:

Remove excess bodily fluid and tissue with a disposable wipe.

Instrument should not be allowed to dry prior to cleaning.

Containment / Transportation:

Follow hospital protocols when handling contaminated and/or biohazardous materials. Instrument should be cleaned within 30 minutes after use to minimize the potential of staining, damage, and drying after use.

Preparation for Decontamination:

If possible, the devices must be reprocessed in a disassembled or opened state.

Pre-Cleaning:

Not required.

Cleaning:

Manual Cleaning Procedure:

1. Prepare a proteolytic enzyme solution according to the manufacturer's instructions.
2. Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with

hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces.

3. Soak instruments for a minimum of 10 minutes. While soaking, scrub surfaces using a nylon-bristled brush until all visible soil has been removed. Actuate moveable mechanisms. Particular attention should be given to crevices, hinged joints, box locks, instrument teeth, rough surfaces, and areas with moving components or springs. Lumens, blind holes, and cannulas should be cleaned using a snug fitting round nylon bristle brush. Insert the snug fitting round brush into the lumen, blind hole, or cannula with a twisting motion while pushing in and out multiple times.

NOTE: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.

4. Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas, and other difficult to access areas.
5. Prepare an ultrasonic cleaning bath with detergent and de-gas according to the manufacturer's recommendations. Completely submerge instruments in the cleaning solution and gently shake them to remove any trapped bubbles. Lumens, blind holes, and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces. Sonically clean the instruments at the time, temperature, and frequency recommended by the equipment manufacturer and optimal for the detergent used. A minimum of ten (10) minutes is recommended.

NOTES:

- **Separate stainless steel instruments from other metal instruments during ultrasonic cleaning to avoid electrolysis.**
 - **Full open hinged instruments**
 - **Use wire mesh baskets or tray designed for ultrasonic cleaners**
 - **Regular monitoring of sonic cleaning performance by means of an ultrasonic activity detector, aluminum foil test, TOSI™ or SonoCheck™ is recommended.**
6. Remove the instruments from the ultrasonic bath and rinse in purified

water for a minimum of one (1) minute or until there is no sign of residue detergent or biologic soil. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas, and other difficult to access areas.

7. Dry Instruments with a clean, absorbent, non-shedding lint free cloth. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas, and difficult to access areas.

Combination Manual/Automated Cleaning Procedure

1. Prepare a proteolytic enzyme solution according to the manufacturer's instructions.
2. Completely submerge instruments in the enzyme solution and gently shake them to remove trapped bubbles. Actuate instruments with hinges or moving parts to ensure contact of the solution with all surfaces. Lumens, blind holes and cannulations should be flushed with a syringe to remove bubbles and ensure contact of the solution with all instrument surfaces.
3. Soak instruments for a minimum of ten (10) minutes. Scrub surfaces using a soft nylon-bristled brush until all visible soil has been removed. Actuate moveable mechanisms. Particular attention should be given to crevices, hinged joints, box locks, instrument teeth, roughened surfaces, and areas with moving components or springs. Lumens, blind holes and cannulas should be cleaned using a snug fitting round nylon bristle brush. Insert the snug fitting round nylon bristle brush into the lumen, blind hole or cannula with a twisting motion while pushing in and out multiple times.

Note: All scrubbing should be performed below the surface of the enzyme solution to minimize the potential of aerosolizing contaminated solution.

4. Remove the instruments from the enzyme solution and rinse in tap water for a minimum of one (1) minute. Actuate all moveable and hinged parts while rinsing. Thoroughly and aggressively flush lumens, holes, cannulas and other difficult to access areas.
5. Place instruments in a suitable washer/disinfector. Follow the washer/disinfector manufacturer's instructions for loading the instruments for maximum cleaning exposure; e.g. open all instruments, place concave instruments on their side or upside down, use baskets and trays designed for washers, place heavier instruments on the bottom of trays and baskets. If the washer/disinfector is equipped with special racks (e.g. for cannulated instruments) use them according to the

manufacturer's instructions.

6. Process instruments using a standard washer/disinfector instrument cycle according to the manufacturer's instructions. The following minimum wash cycle parameters are recommended:

Cycle	Description
1	Pre-wash • Cold Softened Tap Water • 2 minutes
2	Enzyme Spray & Soak • Hot Softened Tap Water • 1 minute
3	Rinse • Cold Softened Tap Water • 1 minute
4	Detergent Wash • Hot Tap Water (64-66°C/146-150°F) • 2 minutes
5	Rinse • Hot Purified Water (64-66°C/146-150°F) • 1 minute
6	Hot Air Dry (116°C/240°F) • 7 – 30 minutes

Notes:

- **The washer/disinfector manufacturer's instructions should be followed.**
- **A washer/disinfector with demonstrated efficacy (e.g. FDA approved, validated to ISO 15883) should be used.**
- **Dry time is shown as a range because it is dependent upon the load size placed into the washer/disinfector.**
- **Many manufacturers pre-program their washer/disinfectors with standard cycles and they may include a thermal low-level disinfection rinse after the detergent wash. The thermal disinfection cycle should be performed to achieve a minimum value A0=600 (e.g. 90°C/194°F for 1 minute according to ISO 15883-1) and is compatible with instruments.**
- **If a lubrication cycle is available that applies to water-soluble lubricant such as Preserve®, Instrument Milk or equivalent material intended for medical device application it is acceptable to use on Instruments unless otherwise indicated.**

Disinfection:

- Instruments and accessories must be terminally sterilized prior to use. See sterilization instructions below.
- Low level disinfection may be used as part of a washer/disinfector cycle but the devices must also be sterilized before use.

Drying:

Dry instruments with a clean, absorbent, non-shedding lint free cloth. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas, and difficult to access areas.

Inspection / Functional Testing:

1. Carefully inspect each device to ensure that all visible blood and soil has been removed.
2. Visually inspect for damage and/or wear.
3. Check the action of moving parts (such as hinges and box-locks) to ensure smooth operation throughout the intended range of motion.
4. Check instruments with long slender features (particularly rotating instruments) for distortion.
5. Where instruments form part of a larger assembly, check that the devices assemble readily with mating components.

***Note: If damage or wear is noted that may compromise the function of the instrument, do not use the instrument, and notify the appropriate responsible person.**

Packaging:

Appropriate packaging for sterilization.

Lubrication and Maintenance:

After cleaning and before sterilization, instruments may be lubricated with a water-soluble lubricant such as Preserve®, Instrument Milk or equivalent material intended for surgical application. Lubricate hinges, threads and other moving part. Always follow the lubricant manufacturer's instructions for dilution, shelf life and application method.

Sterilization:

The Instrument must be cleaned prior to sterilization.

Single devices may be packaged in an approved (e.g. FDA cleared or ISO 11607 compliant) medical grade sterilization pouch or wrap. Care should be used when packaging so that the pouch or wrap is not torn. Devices should be wrapped using the double wrap or equivalent method (ref: AAMI ST79, AORN Guidelines).

Reusable wraps are not recommended.

Moist heat/steam sterilization is the recommended method for the instruments. Use of an approved chemical indicator (class 5) or chemical emulator (class 6) within each sterilization load is recommended. Always

consult and follow the sterilizer manufacturer instructions for load configuration and equipment operation. Sterilizing equipment should have demonstrated efficacy (e.g. FDA clearance, EN 13060 or EN 285 compliance). Additionally the manufacturer's recommendations for installation, validation, and maintenance should be followed.

Validated exposure time and temperature to achieve a 10⁻⁶ sterility assurance level (SAL) is listed below. Local or national specification should be followed where steam sterilization requirements are stricter or more conservative than those listed below.

Method – Pre-Vacuum Steam Sterilization

Cycle Type	Temperature	Exposure Time
United States Recommended Parameters		
Pre-vacuum / Vacuum Pulse	132°C/270°F	4 minutes
Cycle Type	Temperature	Exposure Time
European Recommended Parameters		
Pre-vacuum / Vacuum Pulse	134°C/273°F	3 minutes

Set parameters with four (4) vacuum pulses prior to actual exposure.

A minimum dry time of 60 minutes is recommended. Longer times may be necessary because of load configuration, ambient temperature and humidity, device design, and packaging used.

Flash sterilization is not allowed on lumen instruments!

Storage:

Store sterile, packaged instruments in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes. Note: Inspect every package before use to ensure that the sterile barrier (e.g. wrap, pouch or filter) is not torn, perforated, shows signs of moisture, or appears to be tampered with. If any of those conditions are present then the contents are considered non-sterile and should be reprocessed through cleaning, packaging, and sterilization.

Instructions for Use:

Do not use the instrument if accessories, implant, or other medical devices do not properly fit with the instrument. The correct fit and

geometry of mating accessories, implants or other medical devices will be determined by licensed health care practitioners.

The user of the instrument should ensure that the instrument has no more than 60 minutes of direct exposure to the patient. The user of this set should ensure that no instruments/parts of instruments be left in the patient during and following the surgical procedure.

Inspect device for damage or wear prior to use. Failure to do so could result in the damaging of the device.

***Note: if damage or wear is noted that may compromise the function of the instrument, do not use the instrument and notify the appropriate responsible person.**

The following tables identify the torque ratings of the Male Hex and Male Hexalobe drivers. The drivers should not be subjected to torque loads higher than the identified ratings in order to ensure safe use.

Table 3: Male Hex Driver Torque Rating

Size	Torque Rating (in.-lbs)
1.5 mm	7
2.0 mm	10
2.5 mm	29
7/64"	29
3.0 mm	47
3.2 mm (1/8")	47
3.5 mm	71
5/32"	100
4.5 mm	138
3/16"	138
5.0 mm	224
7/32"	224
1/4"	224
7.0 mm	224

Table 4: Male Hexalobe Driver Torque Ratings

Size	Torque Rating (in.-lbs)
T08	10
T10	17
T15	29
T20	47
T25	71
T27	100
T30	138
T40	240

Disposal:

Devices are reusable; but may be disposed of at the end of the device life by means of standard hospital disposal procedures.

Additional Advices:

It is the duty of the user to ensure that the reprocessing processes are followed; that resources and materials are available to capable personnel; and that hospital protocols are followed. State of the art and often national law requiring these processes and included resources to be validated and maintained properly.

Symbols Used on Labeling¹:

	Caution
	Non-Sterile
	Sterile
	U.S Federal law restricts this device to sale by or on the order of a physician
	CE Mark ¹
	CE Mark with Notified Body # ¹
	Authorized Representative in the European Community
	Manufacturer
	Date of Manufacture
	Lot Number
	Catalog Number
	Consult Instructions for Use
	Medical Device
	Quantity
	Country of Manufacture
Additional Labeling:	
"MANUAL USE ONLY"	Device should not be attached to a power source and is only intended to be handled manually.
"REMOVE FOR CLEANING & STERILIZATION"	Device should be disassembled prior to cleaning and sterilization.

Manufactured By:

Legal Name:
Tecomet, Inc.
5307 95th Avenue
Kenosha, WI 53144 USA

European Representative
Meditec Source GmbH & Co. KG
Sattlerstraße 19 78532
Tuttlingen, Germany

www.tecomet.com

¹Refer to the labelling for CE information

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