



OFFSET CUP IMPACTOR

The Offset Cup Impactor is a reusable hand held instrument designed for use by orthopedic surgeons specifically for the purpose of inserting and implanting acetabular cup components during Total Hip Arthroplasty (THA) surgery.

Modified to fit your custom acetabular cup implant thread connection, our New Product Development Team will evaluate your cup design to engineer the appropriate threaded connection including the required impact sleeve geometry to assure proper fit.

The Offset Cup Impactor features three different ergonomic handle options:

- Black silicone overmolded handle
- Bone Radel® R injection molded handle
- Machined stainless steel handle

The impactor provides a means for holding, inserting, and impacting an acetabular cup implant and includes a custom impact sleeve designed to protect the mating implant interface.

A standard Anteversion Guide (DC087-010-00) with 20 degree anteversion and 45 degrees abduction is included. The Guide attaches to the impactor to provide the surgeon with an external visual guide that replicates the alignment of the instrument to the hip joint during impaction. The anteversion guide is designed to assist during the disassembly process for the cleaning and sterilization process. Please refer to the Tecomet Instruction For Use (IFU) for further instruction and steps required.

Contact your local Sales Manager or visit tecomet.com for more information.



3724 North SR 15
Warsaw, IN 46582 USA
Phone +1 574.267.8700

sales@tecomet.com
www.tecomet.com

Cat #: 5103-00
Rev. E
Printed in the USA ©2023
Tecomet. All rights reserved.