



# PATELLA CALIPER

The **Patella Caliper** is useful for measuring the patella in total knee replacement.

The easy-to-read calibrations range from 0 to 40mm.

Restoring proper patellar thickness is essential for a successful outcome. The Caliper is easily cleaned and features ergonomic handles to provide a secure grip.

Contact your local Sales Manager or visit [tecomet.com](http://tecomet.com) for more information.

Patella Caliper Part #: 14-016-000



3724 North SR 15  
Warsaw, IN 46582 USA  
Phone +1 574.267.8700

[sales@tecomet.com](mailto:sales@tecomet.com)  
[www.tecomet.com](http://www.tecomet.com)

Cat #: 5103-00  
Rev. E  
Printed in the USA ©2023  
Tecomet. All rights reserved.