

STRAIGHT CUP IMPACTOR

The Straight Cup Impactor is a reusable hand held instrument designed for use by orthopedic surgeons specifically for the purpose of manipulating and implanting acetabular cup components during Total Hip Arthroplasty (THA) surgery.

Modified to fit your custom acetabular cup implant thread connection, our New Product Development Team will evaluate your cup design to determine the appropriate fit including the required impact sleeve geometry to assure proper fit. The Straight Cup Impactor features an ergonomic black silicone overmolded handle and includes a Custom Impact Sleeve designed to protect the mating implant interface.

A standard Anteversion Guide (12-047-030) with 15 degree anteversion and 45 degrees abduction is available for purchase. The Guide attaches to the impactor to provide the surgeon with an external visual guide that replicates the alignment of the instrument to the hip joint during impaction.

Contact your local Sales Manager or visit [tecomet.com](https://www.tecomet.com) for more information.

Straight Cup Impactor Part #: 12-051-000



3724 North SR 15
Warsaw, IN 46582 USA
Phone +1 574.267.8700

sales@tecomet.com
www.tecomet.com

Cat #: 5103-00
Rev. E
Printed in the USA ©2023
Tecomet. All rights reserved.