



MODULAR CALCAR PLANER

The **Modular Calcar Planer** provides safe, effective and accurate cutting instrument for the removal of hard calcar bone.

The patented design provides easy disassembly and blade change in a manner of seconds. The Calcar Wrench allows quick turnaround of instrument and intra-operative blade changes.

Surgeons will appreciate a constantly sharp blade producing repeatable, accurate cuts. The disposable blades provide sharp cuts for every procedure.

Our modular blades reduce overall costs by eliminating the need to replace the shaft.

Contact your local Sales Manager or visit tecomet.com for more information.



3724 North SR 15
Warsaw, IN 46582 USA
Phone +1 574.267.8700

sales@tecomet.com
www.tecomet.com

Cat #: 5103-00
Rev. E
Printed in the USA ©2023
Tecomet. All rights reserved.

MODULAR CALCAR PLANER

Part Numbers

*1.5" Calcar Driver Complete Set Options

Complete Set with Hudson Drive End	30-043-000
Complete Set with Large AO Drive End	30-045-000
Complete Set with Zimmer/Tri-Shank Drive End	30-041-000

*1.5" Calcar Driver Complete Set Options

Complete Set with Hudson Drive End	30-044-000
Complete Set with Large AO Drive End	30-046-000
Complete Set with Zimmer/Tri-Shank Drive End	30-042-000

Individual Instrument, Driver Only Options

1.5" Calcar Driver Only, Hudson Drive End	11-103-000
1.5" Calcar Driver Only, Large AO Drive End	11-104-000
1.5" Calcar Driver Only, Zimmer/Tri-Shank Drive End	11-074-000
1.8" Calcar Driver Only, Hudson Drive End	11-105-000
1.8" Calcar Driver Only, Large AO Drive End	11-106-000
1.8" Calcar Driver Only, Zimmer/Tri-Shank Drive End	11-076-000

Disposable Calcar Blades

1.5" Disposable Calcar Blade	11-075-000
1.8" Disposable Calcar Blade	11-077-000

Calcar Driver Wrench

Calcar Driver T-Handle Wrench	11-078-000
-------------------------------	------------

*Complete sets include: one modular calcar planer, one disposable calcar blade and one calcar driver t-handle wrench.

